

FIWA®

product catalogue

our product portfolio for hospitals
either in the areas of ward and surgery,
for outpatient facilities and
for medical practices

COMPRESSES

FIWA®komp	• Gauze swabs, non-sterile	9
	• Gauze swabs, sterile	11
	• Gauze swabs, sterile (double-packed)	13
	• Gauze swabs X-RAY, non-sterile	15
	• Gauze swabs X-RAY, sterile	17
	• Swabs Y-cut GAUZE, sterile	19
	• Ear swabs, sterile	21
	• Rolled swabs X-RAY, sterile	23
FIWA®vlies	• Nonwoven swabs, non-sterile	25
	• Nonwoven swabs, sterile	27
	• Swabs Y-cut NONWOVEN, non-sterile	29
	• Swabs Y-cut NONWOVEN, sterile	31
	• Cellulose-nonwoven-swabs, non-sterile	33
FIWA®sorb	• Absorbent pads GAUZE, non-sterile	35
	• Absorbent pads GAUZE, sterile	37
	• Absorbent pads NONWOVEN, non-sterile	39
	• Absorbent pads NONWOVEN, sterile	41

DRESSING GAUZE

FIWA®mull	• Dressing gauze, non-sterile	43
------------------	-------------------------------------	----

GAUZE BALLS

FIWA®sling	• Gauze balls, non-sterile	46
	• Gauze balls, sterile	48
	• Gauze balls X-RAY, non-sterile	50
	• Gauze balls X-RAY, sterile	52
	• Gauze balls, security thread, non-sterile	54
	• Gauze balls, security thread, sterile	56
	• Pointed gauze balls X-RAY, non-sterile	58
	• Pointed gauze balls X-RAY, sterile	60
	• Pointed gauze balls, security thread X-RAY, sterile	62
FIWA®prep	• Preparation balls X-RAY, non-sterile	64
	• Preparation balls X-RAY, sterile	66
FIWA®zell	• Cellulose squares, non-sterile	68
	• Cellulose squares germ-reduced, non-sterile	70
	• Dispenser for cellulose swabs, non-sterile	72

LAPORATOMY SPONGES

FIWA®lap	• Laparotomy sponges X-RAY, non-sterile	73
	• Laparotomy sponges X-RAY, sterile	76
	• Laparotomy sponges nonwoven X-RAY, sterile	79

INTESTINAL BAGS

FIWA®intest	• Intestinal bag X-RAY, non-sterile	81
	• Intestinal bag X-RAY, sterile	83

TAMPONDADES

FIWA®tamp	• Gauze tamponade, sterile	85
	• Gauze tamponade X-RAY, sterile	86

MEDICAL APPLICATORS

FIWA®app	• Medical applicators, non-sterile	88
	• Medical applicators, sterile	90

BANDAGES (NOT ADHESIVE)

FIWA®elast	• Conforming bandages, non-sterile	92
	• Conforming bandages, sterile, (single packed)	95
	• Conforming bandages, sterile (double-packed)	97
FIWA®elast krepp	• Conforming bandages CREPE, non-sterile	99
FIWA®in-elast	• Gauze bandages, non-sterile	101
FIWA®uni	• Elastic universal-bandage, non-sterile	103
	• Elastic universal-bandage, sterile, (double-packed)	105
FIWA®id-uni	• Ideal universal bandages, non-sterile	106
FIWA®ideal	• Ideal bandages, non-sterile	108
	• Ideal bandages, sterile, (single packed)	110
	• Ideal bandages, sterile, (double-packed)	112

	• Ideal bandages coloured, non-sterile	114
FIWA®BW-plus	• Ideal bandages BW-plus, non-sterile	116
FIWA®ideal DIN	• Ideal bandages DIN 61632, non-sterile	118
	• Ideal bandages DIN 61632, sterile, (single packed)	120
	• Ideal bandages DIN 61632, sterile, (double-packed)	122
FIWA®dur fein	• Long-stretch bandages, non-sterile	124
FIWA®lan	• Short-stretch bandages, non-sterile	126
	• Short-stretch bandages, sterile, (double-packed)	128
FIWA®lan twin	• Short-stretch bandages, non-sterile	129
FIWA®lan ultra	• Ultra-short-stretch bandages, non-sterile	131
 COHESIVE BANDAGES		
FIWA®flex haft	• Cohesive dressing bandage, non-sterile	133
FIWA®elast haft	• Cohesive conforming bandages, non-sterile	135
	• Cohesive conforming bandages coloured, non-sterile	137
FIWA®elast krepp haft	• Cohesive conforming bandages CREPE, non-sterile	139
	• Cohesive conforming bandages CREPE coloured, non-sterile	141
FIWA®ideal haft	• Cohesive ideal bandages, non-sterile	143
FIWA®lan haft	• Cohesive short-stretch bandages, non-sterile	145
 ADHESIVE BANDAGES		
FIWA®acryl	• Acrylic bandage, non-sterile	147
FIWA®bi-elast	• Adhesive bandage, non-sterile	149

SPECIAL BANDAGES

FIWA®soft	• Padding bandages, non-sterile	151
FIWA®crepex	• Crepex bandages, sterile	153
FIWA®krüll	• Krinkle gauze rolls, non-sterile	155
	• Krinkle gauze rolls, sterile	157
FIWA®zink	• Zinc paste bandage, non-sterile	158

BANDAGE CLIPS

FIWA®clip	• Bandage clips, non-sterile	160
------------------	------------------------------------	-----

DRESSING

FIWA®net	• Tubular net dressing, non-sterile	161
FIWA®tube	• Tube dressing, non-sterile	163
	• Tube finger-dressing, non-sterile	165
FIWA®tube fast	• Tube dressing, longitudinal and transverse elastic, non-sterile	167
FIWA®tape	• Sports tape, non-sterile	169
FIWA®kin	• Elastic tape bandage for kinesiological therapy, non-sterile	171

PLASTER PRODUCTS

FIWA®med steril	• Adhesive wound dressing, nonwoven, sterile	173
	• Adhesive wound dressing, transparent, sterile	175
FIWA®strip	• Injection strips, non-sterile	177
	• Finger plaster, non-sterile	179
FIWA®stretch	• Nonwoven retention tape, non-sterile	181

Table of contents

FIWA®stretch foil	• Foil dressing, transparent, non-sterile	183
FIWA®med soft	• Woundplaster NONWOVEN, non-sterile	185
FIWA®med tex	• Woundplaster TEXTILE, non-sterile	187
FIWA®silk	• Medical tape SILK, non-sterile	189
FIWA®tex	• Medical tape TEXTILE, non-sterile	191
FIWA®por	• Medical tape PORE, non-sterile	193
FIWA®clear	• Medical tape TRANSPARENT, non-sterile	195

CELLULOSE/HYGIENIC PRODUCTS

FIWA®protect	• Cellulose, non-sterile	197
	• Underpads, non-sterile	209
	• Protective covers, non-sterile	201
	• Protective covers reinforced, non-sterile	203
	• Dialysis underpad, non-sterile	205
	• Multi-purpose bag, non-sterile	207

PROTECTION PRODUCTS

FIWA®protect	• Surgical facemasks, non-sterile	208
---------------------	---	-----

SET-SYSTEMS

FIWA®set	• Sets, sterile	210
	• Sets X-RAY, sterile	212

MISCELLANEOUS

Dr. Till	• Bag	214
kids emergency box	• Refill kit	216
finellis	• Mini sanitary bag for tampons	218

Product description:

17-thread swabs made of gauze according to DIN EN 14079,
100 % cotton,
with tucked-in cut edges

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, to clean the skin and wound surface, as dressing, for haemostasis, only for external use

Product image (example):



Product variants:

	SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
8-ply	5 x 5 cm	731015	16657529	4260604533632	10	variable
	5 x 5 cm	731011	08821400	4260348610088	100	10
	7,5 x 7,5 cm	731035	16657535	4260604533649	10	variable
	7,5 x 7,5 cm	731031	08821423	4260348610095	100	10
	10 x 10 cm	731065	16657541	4260604533656	10	variable
	10 x 10 cm	731061	08821452	4260348610101	100	10
	10 x 20 cm	731091	08821481	4260348610118	100	10
12-ply	5 x 5 cm	731021	08821417	4260348610125	100	10
	7,5 x 7,5 cm	731041	08821446	4260348610132	100	10
	10 x 10 cm	731071	08821469	4260348610149	100	10
	10 x 20 cm	731101	08821498	4260348610156	100	10

Product data sheet

	SIZE	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
16-ply	5 x 5 cm	731024	04808319	4260348610163	100	10
	7,5 x 7,5 cm	731051	04797254	4260348610170	100	10
	10 x 10 cm	731081	08821475	4260348610187	100	10
24-ply	10 x 20 cm	731104	-	4260604532802	50	10

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
 The product has to be sterilised using a validated procedure before coming into contact with open wounds.
 The product has not to be used, if the packaging is damaged.
 The product has not to be reused.
 The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)
 Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)
 The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-01

Product description:

17-thread swabs made of gauze according to DIN EN 14079,
100 % cotton,
with tucked-in cut edges

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, to clean the skin and wound surface, as dressing, for haemostasis, only for external use

Product image (example):



Product variants:

STERILE sealed (single packaged):

	SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
8-ply	5 x 5 cm	321015	04676817	4260348610194	5 x 2	variable
		321013	02773775	4260348610200	25 x 2	40
		321014	02135738	4260348610217	20 x 5	20
	7,5 x 7,5 cm	321035	04676823	4260348610224	5 x 2	variable
		321033	02773798	4260348610231	25 x 2	30
	10 x 10 cm	321069	04676846	4260348610248	5 x 2	variable
		321063	02773806	4260348610255	25 x 2	20
		321068	02135885	4260348610262	20 x 5	12
		321065	10822720	4260348610279	10 x 10	12
	10 x 20 cm	321094	04676852	4260348610286	5 x 2	variable
		321092	04317998	4260348610293	25 x 2	12

Product data sheet

	SIZE	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
12-ply	5 x 5 cm	321023	02137192	4260348610309	25 x 2	60
		321024	11859872	4260348610316	20 x 5	16
	7,5 x 7,5 cm	321043	02137200	4260348613850	25 x 2	30
		321044	05919073	4260348610354	20 x 5	30
	10 x 10 cm	321076	02141331	4260348610361	25 x 2	12
		321074	02143778	4260348610378	20 x 5	12
		321075	02143979	4260348610385	10 x 10	12
	10 x 20 cm	321103	02141348	4260348610408	25 x 2	12
		321101	10822743	4260348610415	20 x 5	12
321105		10822766	4260348610422	10 x 10	12	
16-ply	10 x 10 cm	321083	02144105	4260348610439	25 x 2	12
		321084	-	4260348610446	20 x 5	12
		321085	02144140	4260348610453	10 x 10	12

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Is / 4 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-01

Product description:

17-thread swabs made of gauze according to DIN EN 14079,
100 % cotton,
with tucked-in cut edges

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, to clean the skin and wound surface, as dressing, for haemostasis

Product image (example):



Product variants:

STERILE sealed (double packed with counting cards for the operating room):

	SIZE	PCS/ POUCH	X-RAY	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
8-ply	10 x 10 cm (17-thread)	25	NO	520071	-	4260604533588	24	4
12-ply	10 x 10 cm (17-thread)	10	NO	520030*	-	4260348614581	30	4
16-ply	10 x 10 cm (17-thread)	10	NO	520001	-	4260348614666	40	4
32-ply	10 x 10 cm (17-thread)	10	NO	520060*	-	4260604533762	18	4

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Is / 4 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

*Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-01

Product description:

17-thread swabs made of gauze according to DIN EN 14079,
100 % cotton,
with tucked-in cut edges,
with X-RAY contrast thread

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, to clean the skin and wound surface, as dressing, for haemostasis

Product image (example):



Product variants:

	SIZE	X-RAY	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
8-ply	10 x 10 cm	with	731062	-	4260604532277	100	10
12-ply	5 x 5 cm	with	731023	11860651	4260348614444	100	10
	7,5 x 7,5 cm	with	731042	11860668	4260348614451	100	10
	10 x 10 cm	with	731072	11860680	4260348614468	100	10
	10 x 10 cm GREEN <i>bundled à 5 pieces</i>	with	731074	-	4260348614475	100	40
	10 x 13 cm	with	731122	-	4260348614482	100	10
	10 x 20 cm	with	731102	-	4260348614499	100	10

Product data sheet

	SIZE	X-RAY	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
16-ply	5 x 5 cm	with	731022	11860645	4260348614505	100	10
	7,5 x 7,5 cm	with	731052	11860674	4260348614512	100	10
	10 x 10 cm	with	731082	02920058	4260348614529	100	10
	10 x 20 cm	with	731114	-	4260348614536	100	10
32-ply	10 x 10 cm	with	731112	15204990	4260348614543	100	10

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
 The product has to be sterilised using a validated procedure before coming into contact with open wounds.
 The product has not to be used, if the packaging is damaged.
 The product has not to be reused.
 The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: IIa / 7 (annex IX directive 93/42/EEC)
 Conformity assessment: annex VII + V (directive 93/42/EWG)
 The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-01

Product description:

17-, 20- or 24-thread swabs made of gauze according to DIN EN 14079,
100 % cotton,
with tucked-in cut edges,
with X-RAY contrast thread

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, to clean the skin and wound surface, as dressing, for haemostasis

Product image (example):



Product variants:

STERILE sealed (double packed with counting cards for the operating room):

	SIZE	PCS/ POUCH	X-RAY	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
12-ply	7,5 x 7,5 cm (17-thread)	10	YES	520007	11518214	1234567890123	80	4
		20	YES	520008	-	4260348614567	45	4
	10 x 10 cm (17-thread)	10	YES	520009	11860384	4260348614574	60	4
		20	YES	520010	11518220	4260348614598	36	4
		50	YES	520011	-	4260348614604	10	4
	10 x 13 cm (17-thread)	40	YES	520013	-	4260348614611	12	4
	10 x 20 cm (17-thread)	10	YES	520015	11860409	4260348614628	40	4
		20	YES	520016	-	4260348614635	20	4

Product data sheet

	SIZE	PCS/ POUCH	X-RAY	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
16-ply	5 x 5 cm (17-thread)	20	YES	520006	-	4260348614642	60	4
	10 x 10 cm (17-thread)	10	YES	520002	11860361	4260348614659	40	4
		20	YES	520003	11860378	4260348614673	30	4
		40	YES	520004	-	4260348614680	15	4
	10 x 13 cm (20-thread)	20	YES	520014	-	4260348614697	16	4
	10 x 20 cm (17-thread)	10	YES	520038	-	4260348614703	30	4
32-ply	10 x 10 cm (17-thread)	10	YES	520012	11860390	4260348614710	30	4
		10	YES (over long page)	520061	-	4260348616394	24	4
Edge-stitched:								
8-ply	12 x 38 cm (20-thread)	10	YES	520062	-	4260348614741	12	4
GREEN, edge-stitched, pre-washed:								
2-ply	10 x 10 cm (24-thread)	10	YES	231310	-	4260348614758	30	4

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: IIa / 7 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-01

Product description:

17-thread swabs with a Y-shaped cut,
made of gauze according to DIN EN 14079,
100 % cotton,
with tucked-in cut edges

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, as dressing, for haemostasis, only for external use

Product image (example):



Product variants:

STERILE sealed (single packaged):

	SIZE	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
8-ply	5 x 5 cm	324033	01023345	4260348610477	25 x 2	40
	7,5 x 7,5 cm	324013	01023339	4260348610484	25 x 2	30
	10 x 10 cm	324023	01023322	4260348610491	25 x 2	20
12-ply	5 x 5 cm	324034	15436424	4260348610507	25 x 2	40
	7,5 x 7,5 cm	324014	11859984	4260348610514	25 x 2	30
	10 x 10 cm	324024	15436418	4260348610521	25 x 2	12

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Is / 4 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-04

Product description:

17-thread swabs with a semicircular cut-out for better application behind the ear made of gauze according to DIN EN 14079, 100 % cotton, with tucked-in cut edges

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, as dressing, for haemostasis, only for external use

Product image (example):

Product variants:

STERILE sealed (single packaged):

	SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
32-ply	5 x 10 cm	315012	-	4260604532857	25 x 2	12

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Is / 4 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-01

Product data sheet

Product description:

rolled swabs made of gauze according to DIN EN 14079,
 100 % cotton,
 with tucked-in cut edges,
 wrapped with thread,
 with X-RAY contrast thread

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, for haemostasis, particularly suitable for minimally invasive operations with a Da Vinci surgical system

Product image (example):



Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE	X-RAY	REF	PZN (POUCH)	GTIN	PCS/POUCH	POUCH/BOX
10 cm of length 1,2 cm diameter	YES	520068	14130545	4260348616769	5	50
5 cm of length 1,5 cm diameter	YES	520073	-	4260604539481	5	20

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-01

Product description:

non-woven swabs (30 g/m²) made of 70 % rayon and 30 % polyester fibres,
open, gauze-like structure

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, to clean the skin and wound surface, as dressing, for haemostasis, only for external use

Product image (example):



Product variants:

	SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
4-ply	5 x 5 cm	1727815001	02249241	4260348610538	200	20
	7,5 x 7,5 cm	1727815002	02249235	4260348610545	200	20
	10 x 10 cm	211091	02249229	4260348610552	200	20
	10 x 20 cm	1727815004	04808911	4260348610569	200	10
6-ply	5 x 5 cm	221091	04808041	4260348610576	100	10
	7,5 x 7,5 cm	221092	04808058	4260348610583	100	10
	10 x 10 cm	221093	04808070	4260348610590	100	10
	10 x 20 cm	221094	04808207	4260348610606	100	10

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-03

Product description:

non-woven swabs (30 g/m²) made of 70 % rayon and 30 % polyester fibres,
open, gauze-like structure

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, to clean the skin and wound surface, as dressing, for haemostasis, only for external use

Product image (example):



Product variants:

STERILE sealed (single packaged):

	SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
4-ply	5 x 5 cm	328011	04808271	4260348610613	25 x 2	100
	7,5 x 7,5 cm	328021	04808302	4260348610620	25 x 2	50
	10 x 10 cm	328031	04793807	4260348610637	25 x 2	40
	10 x 20 cm	328041	04793813	4260348610644	25 x 2	10
6-ply	5 x 5 cm	328111	04793836	4260348610651	25 x 2	100
	7,5 x 7,5 cm	328121	04793842	4260348610675	25 x 2	50
		328122	11860332	4260348610682	20 x 5	40
	10 x 10 cm	328131	04793859	4260348610699	25 x 2	30
		328132	11860349	4260348610705	20 x 5	16
		328133	11860355	4260348610712	24 x 25	4
10 x 20 cm	328141	04793865	4260348610729	25 x 2	10	

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Is / 4 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-03

Product description:

non-woven swabs (30 g/m²) with a Y-shaped cut made of 70 % rayon and 30 % polyester fibres, open, gauze-like structure

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, as dressing, for haemostasis, only for external use

Product image (example):

Product variants:

	SIZE	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
4-ply	7,5 x 7,5 cm	221096	11859731	4260348610736	200	40

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-03

Product description:

non-woven swabs (30 g/m²) with a Y-shaped cut made of 70 % rayon and 30 % polyester fibres, open, gauze-like structure

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, as dressing, for haemostasis, only for external use

Product image (example):



Product variants:

STERILE sealed (single packaged):

	SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
4-ply	7,5 x 7,5 cm	324114	11860036	4260348610743	25 x 2	50
	10 x 10 cm	324125	11860042	4260348610750	25 x 2	40
6-ply	5 x 5 cm	324102	07200529	4260348610767	25 x 2	100
	7,5 x 7,5 cm	324112	07200506	4260348610781	25 x 2	50
	10 x 10 cm	324123	07200340	4260348610798	25 x 2	30

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Is / 4 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-03

Product description:

wound side: non-woven fabric made of 100 % polypropylene,
absorbent layer: highly bleached cellulose (according to DAB),
both layers connected by two seams

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, to clean the skin and wound surface, as dressing, for haemostasis, only for external use

Product image (example):



Product variants:

SIZE	REF	PZN (BAG)	GTIN	ROLLS/ BAG	BAGS/ CARTON
10 cm x 10 m	1739023	01902443	4260348613942	1	18

Warnings / contraindications:

- The product is not to be used in the case of an existing allergy to the material.
- The product has to be sterilised using a validated procedure before coming into contact with open wounds.
- The product has not to be used, if the packaging is damaged.
- The product has not to be reused.
- The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-03

Product description:

absorbent cotton wool core made of a cotton and rayon blend,
soft, natural and skin-friendly gauze cover made of 17-thread gauze according to DIN EN 14079
(100 % cotton), reverse side made of a water-repellent but air-permeable layer

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, as dressing, for haemostasis, only for external use

Product image (example):

Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
15 x 20 cm	432042	-	4260604539061	20	variable
20 x 30 cm	432041	07200280	4260348610804	20	20

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-05


Product description:

absorbent cotton wool core made of a cotton and rayon blend,
soft, natural and skin-friendly gauze cover made of 17-thread gauze
according to DIN EN 14079 (100 % cotton),
reverse side made of a water-repellent but air-permeable layer

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, as dressing, for haemostasis, only for external use

Product image (example):



Product variants:

STERILE sealed (single packaged):

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
20 x 30 cm	322041	06727743	4260348610811	12 x 1	9

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Is / 4 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-39-04-00

Product description:

2 centre layers made of liquid-spreading tissue fabric,
 core made of cellulose flakes,
 soft, skin-friendly fleece cover,
 waterproof reverse side made of a blue polypropylene layer

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, as dressing,
 for haemostasis, only for external use

Product image (example):

Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
10 x 10 cm	132020	07200311	4260348610828	25	63
10 x 20 cm	132021	07200305	4260348610835	25	16
20 x 20 cm	132031	07200297	4260348610842	25	8
20 x 40 cm	132041	06864003	4260348610859	25	8

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-05

Product description:

2 centre layers made of liquid-spreading tissue fabric,
core made of cellulose flakes,
soft, skin-friendly fleece cover,
waterproof reverse side made of a blue polypropylene layer

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, as dressing, for haemostasis, only for external use

Product image (example):



Product variants:

STERILE sealed (single packaged):

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
10 x 10 cm	322014	14367282	4260604532352	10 x 1	54
	322011	00527658	4260348610866	50 x 1	9
10 x 20 cm	322024	14367299	4260604532413	10 x 1	48
	322021	00527670	4260348610873	25 x 1	9
15 x 25 cm	322033	14411586	4260348613973	12 x 1	9
20 x 20 cm	322031	00527724	4260348610880	12 x 1	9
20 x 30 cm	322043	11859926	4260348610897	12 x 1	9
20 x 40 cm	322032	06863989	4260348610903	10 x 1	9

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Is / 4 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-02-05

Product description:

17- or 20-thread gauze according to DIN EN 14079,
100 % cotton

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, to clean the skin and wound surface, as dressing, for haemostasis, only for external use

Product image (example):



Product variants:

rolled (20-threads): 80 cm wide layer laid to 10 cm width ⇒ 8-ply

SIZE	REF	PZN (BOX)	GTIN	ROLLS/ BOX	BOXES/ CARTON
1m x 10cm	1708815001	08821788	4260348611139	1	100
2m x 10cm	1708815002	08821794	4260348611146	1	100
5m x 10cm	1708815003	08821802	4260348611153	1	20
10m x 10cm	1708815004	08821819	4260348611160	1	20
40m x 10cm	1708815009	08821825	4260348611177	1	10



zigzag (20-threads): 80 cm wide layer laid to 10 cm width ⇒ 8-ply

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
1m x 10cm	1708815005	08821765	4260348611184	1	120
2m x 10cm	1708815006	08821759	4260348611191	1	80
5m x 10cm	1708815007	08821742	4260348611207	1	40
10m x 10cm	1708815008	08821736	4260348611214	1	20

cut (20-threads):

SIZE	REF	PZN (BAG)	GTIN	PCS/ BAG	BAGS/ CARTON
25 x 25 cm	1704815001	11859665	4260348611238	500	18
27 x 32 cm	250515	-	4260348616936	500	10

bales (17-threads):

SIZE	REF	PZN (BAG)	GTIN	METER/ BAG	BAGS/ CARTON
64 cm x 100 m	1406417	10763621	4260348611245	100	10

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-09-04

Product description:

20-thread ready-to-use gauze balls made of gauze according to DIN EN 14079,
100 % cotton

Field of application / indication:

for care of secreting wounds, for sucking up blood or other body fluids, to clean the skin and wound surface,
for haemostasis, only for external use

Product image (example):



Product variants:

in a PE-Bag:

SIZE	REF	PZN (BAG)	GTIN	PCS/BAG	BAGS/CARTON
8 x 8 cm, hazelnut-sized	733021	08821601	4260348610989	500	20
15 x 15 cm, walnut-sized	733041	08821618	4260348610996	500	16
20 x 20 cm, plum-sized	733051	08821624	4260348611009	500	10
25 x 25 cm, egg-sized	733061	08821630	4260348611016	500	8
30 x 30 cm, (extra-)large	733071	08821647	4260348611023	500	6
30 x 40 cm, extra-large	733081	08821653	4260348611030	500	6

in a quick-release bag with holes (for hanging):

SIZE	REF	PZN (BAG)	GTIN	PCS/BAG	BAGS/CARTON
20 x 20 cm, plum-sized	733054	-	4260348616820	50	100

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
The product has to be sterilised using a validated procedure before coming into contact with open wounds.
The product has not to be used, if the packaging is damaged.
The product has not to be reused.
The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)
Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)
The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-01-01

Product description:

20-thread ready-to-use gauze balls made of gauze according to DIN EN 14079,
 100 % cotton

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, to clean the skin and wound surface, for haemostasis, only for external use

Product image (example):

Product variants:

STERILE sealed (single packaged):

* The gauze balls can be removed from both sides of the bag without having to tear open the rest of the bag.
 Of course, it can also be torn open completely and all the gauze balls removed at once.

SIZE	PCS/ POUCH	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
20 x 20 cm plum-sized	2	323103	16657587	4260604533670	10 x 2	variable
	2	238025	08821682	4260348611047	50 x 2	12
	2+2*	323097	14305772	4260604532383	25 x 4	12
	2+3*	323055	04767276	4260348611054	30 x 5	12
	10	323089	-	4260348616899	20 x 10	12
	25	323056	04767299	4260348611061	10 x 25	12

Product data sheet

25 x 25 cm egg-sized	4	323071	11859961	4260348611078	35 x 4	12
	5	323074	11859978	4260348611085	25 x 5	12
	2+3*	323057	00406133	4260348611092	25 x 5	12
30 x 40 cm extra-large	5	323059	11859955	4260348611108	30 x 5	4

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Is / 4 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-01-01

Product description:

20-thread ready-to-use gauze balls made of gauze according to DIN EN 14079,
100 % cotton,
with X-RAY contrast thread

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, to clean the skin and wound surface, for haemostasis

Product image (example):



Product variants:

packed in PE-bags:

SIZE	X-RAY	REF	PZN (BAG)	GTIN	PCS/ BAG	BAGS/ CARTON
11 x 11 cm	YES	733032	-	4260348614765	500	16
15 x 15 cm	YES	733042	-	4260348614772	500	16
20 x 20 cm	YES	733052	11860705	4260348614789	500	10
25 x 25 cm	YES	733062	11860711	4260348614796	500	8
30 x 30 cm	YES	733072	-	4260348614802	500	4
30 x 40 cm	YES	733082	11860728	4260348614819	500	6
40 x 50 cm	YES	740718	11860734	4260348614826	500	4

sliding box (SB) with grid:

SIZE	X-RAY	REF	PZN (BOX)	GTIN	PCS/SB	SB/BOX	BOXES/CARTON
25 x 25 cm	YES	733067	-	4260604539191	10	50	6
25 x 25 cm	YES	733068	-	4260604539085	5	50	6
30 x 40 cm	YES	733069	-	4260604539108	10	50	3

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: IIa / 7 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-01-01

Product description:

20-thread ready-to-use gauze balls made of gauze according to DIN EN 14079,
100 % cotton,
with X-RAY contrast thread

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, to clean the skin and wound surface, for haemostasis

Product image (example):



Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE	PCS/ POUCH	X-RAY	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
11 x 11 cm	10	YES	560057*	-	4260348614833	35	8
20 x 20 cm	10	NO	560008	-	4260348614840	60	4
	10	YES	560001	11533633	4260348614857	60	4
25 x 25 cm	10	YES	560002	11860496	4260348614864	45	4
	20	YES	560006	-	4260348614871	30	4
30 x 30 cm	10	YES	560003	-	4260348614888	30	4
	20	YES	560007	-	4260348614895	28	4
30 x 40 cm	10	YES	560004	11860504	4260348614901	30	4
	20	YES	560016	-	4260348614918	15	4
40 x 50 cm	10	YES	560005	11860510	4260348614925	15	4

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: IIa / 7 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

*Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-01-01

Product description:

20-thread ready-to-use gauze balls made of gauze according to DIN EN 14079,
100 % cotton,
with a 30 cm long security thread

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, for haemostasis,
only for external use

Product image (example):

Product variants:

SIZE	REF	PZN (BAG)	GTIN	PCS/ BAG	BAGS/ CARTON
15 x 15 cm, walnut-sized	250139	7200334	4260348614932	100	variable
25 x 25 cm, egg-size	250163	7200328	4260348614949	100	variable

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

The security thread is used to locate the swab and not to remove it. The product may only be removed from the surgical field using surgical instruments.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-01-01

Product description:

20-thread ready-to-use gauze balls made of gauze according to DIN EN 14079,
 100 % cotton,
 with a 30 cm long security thread,
 with X-RAY contrast thread

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, for haemostasis,
 only for external use

Product image (example):

Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE	PCS/ POUCH	X-RAY	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
15 x 15 cm, walnut-sized	5	YES	560055	11860533	4260348614956	30	5
20 x 20 cm, plum-sized	10	YES	560039	11860527	4260348614963	25	6
25 x 25 cm, egg-sized	5	YES	560056	11860562	4260348614970	25	5

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

The security thread is used to locate the product and not to remove it. The product may only be removed from the surgical field using surgical instruments.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-01-01


Product description:

20-thread ready-to-use pointed gauze balls made of gauze according to DIN EN 14079,
100 % cotton,
conical shape,
made from one piece of gauze,
with X-RAY contrast thread

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, for haemostasis

Product image (example):



Product variants:

SIZE	X-RAY	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
6 x 8 cm	YES	250211	08821860	4260348615199	500	1
15 x 15 cm	YES	250214	08821883	4260348615205	500	1

Warnings / contraindications:

- The product is not to be used in the case of an existing allergy to the material.
- The product has to be sterilised using a validated procedure before coming into contact with open wounds.
- The product has not to be used, if the packaging is damaged.
- The product has not to be reused.
- The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: IIa / 7 (annex IX directive 93/42/EEC)

Conformity assessment: annex II without section (4) (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-01-03

Product description:

20-thread ready-to-use pointed gauze balls made of gauze according to DIN EN 14079,
 100 % cotton,
 conical shape,
 made from one piece of gauze,
 with X-RAY contrast thread

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, for haemostasis

Product image (example):



Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE	PCS/ POUCH	X-RAY	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
6 x 8 cm	10	YES	580008	-	4260348615212	50	variable
15 x 15 cm	6	YES	580009	11860616	4260348615243	30	variable

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-01-03

Product description:

20-thread ready-to-use pointed gauze balls made of gauze according to DIN EN 14079,
100 % cotton,
conical shape,
made from one piece of gauze,
with a 30 cm long security thread

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, for haemostasis

Product image (example):



Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE	PCS/ POUCH	X-RAY	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
6 x 8 cm	10	YES	580015	-	4260348615229	50	variable
15 x 15 cm	6	YES	580003	-	4260348616905	30	12

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

The security thread is used to locate the pointed swab and is not suitable for its removal. The pointed swabs may only be removed from the surgical field using surgical instruments.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-01-03

Product description:

24-thread ready-to-use preparation balls made of gauze according to DIN EN 14079,
 100 % cotton,
 firmly twisted,
 made from one piece of gauze,
 with X-RAY contrast thread

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, for haemostasis

Product image (example):

Product variants:

packed in PE-bags:

SIZE	X-RAY	REF	PZN (BAG)	GTIN	PCS/BAG	BAGS/CARTON
6 x 6 cm	YES	733111	08821831	4260348614987	1000	variable
8 x 8 cm	YES	733121	08821848	4260348614994	1000	variable
12 x 12 cm	YES	733131	08821854	4260348615007	1000	variable

empty sliding boxes (SB):

SIZE	PCS/SB	REF	PZN	GTIN	PCS/SB	SB/CARTON
5 x 9 cm	e.g. for 10 preparation balls	733150	-	4260348615014	variable	variable

sliding box (SB) with grid:

SIZE	PCS/ SB	X-RAY	REF	PZN	GTIN	SB/ BOX	BOXES/ CARTON
6 x 6 cm	10	YES	733112	-	4260348615021	100	20
	20	YES	733160	-	4260348615038	100	20
8 x 8 cm	10	YES	733116	-	4260348615045	100	20
	20	YES	733157	-	4260348615052	100	20
10 x 10 cm	10	YES	733118	-	4260348615069	100	20
	20	YES	733161	-	4260348615076	50	10
12 x 12 cm	10	YES	733120	-	4260348615083	50	10
	20	YES	733159	-	4260348615090	20	10

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: IIa / 7 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-01-02

Product description:

24-thread ready-to-use preparation balls made of gauze according to DIN EN 14079,
100 % cotton,
firmly twisted,
made from one piece of gauze,
with X-RAY contrast thread

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, for haemostasis

Product image (example):



Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE	PCS/ POUCH	X-RAY	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
6 x 6 cm	10	YES	570001	11860579	4260348615106	60	4
6 x 6 cm	20	YES	570004	-	4260348615113	60	4
8 x 8 cm	10	YES	570002	11860585	4260348615120	60	4
8 x 8 cm	20	YES	570010	-	4260348615137	60	4
8 x 8 cm	30	YES	570007	-	4260348615144	60	4
10 x 10 cm	10	YES	570009	-	4260348615151	60	4
12 x 12 cm	10	YES	570003	11860591	4260348615168	60	4
12 x 12 cm	20	YES	570006	-	4260348615175	42	4

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: IIa / 7 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-01-02

Product data sheet

Product description:

roll with cellulose swabs arranged in pairs made of highly bleached cellulose according to DAB,
220 g/m²,
perforated for tearing off

Field of application / indication:

for sucking up blood or other body fluids, to clean the skin, only for external use

Product image (example):



Product variants:

SIZE	REF	PZN (BAG)	GTIN	PCS/ ROLLS	ROLLS/ BAG	BAGS/ CARTON
4 x 5 cm, 2 rolls à 500 pcs.	1739021	08822718	4260348611115	500	2	16

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-22-12-13

Product description:

roll with cellulose swabs arranged in pairs made of highly bleached cellulose according to DAB,
220 g/m²,
perforated for tearing off,
germ-reduced (treated with ethylene oxide)

Field of application / indication:

for sucking up blood or other body fluids, to clean the skin, only for external use

Product image (example):

Product variants:

SIZE	REF	PZN (BAG)	GTIN	PCS/ ROLL	ROLLS/ BAG	BAGS/ CARTON
4 x 5 cm, 1 roll à 500 pcs.	1739025	05103288	4260348611122	500	1	32

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product is non-sterile but germ-reduced. It was treated with ethylene oxide.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-22-12-13

Product description:

transparent dispenser with rubber feet

Field of application / indication:

Suitable as a dispenser for the hygienic storage of rolled products, such as cellulose swabs, cellulose non-woven compresses or rolled gauze

Product variants:

DIAMETER	WIDTH	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
approx. 20 cm	approx. 10,5 cm	1739022	-	4260348616370	1	variable

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

The product does not contain any substances according to REACH.

eCI@ss:

34-21-01-92

Product data sheet

Product description:

20- or 28-thread gauze cloth according to DIN EN 14079, 100 % cotton, with or without loop (30 cm), pre-washed (unless otherwise stated), with X-RAY contrast chip

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, for haemostasis

Product image (example):



Product variants:

SIZE	PLY	COLOUR	X-RAY	LOOP	REF	PZN (BOX) GTIN	PCS/ BUNDLE	BUNDLE/ BOX	BOXES/ CARTON
8 x 90 cm	4-ply	green	YES	YES	29196	11859843 4260348615267	4 x 5	8	4
				NO	29194	- 4260348615274	4 x 5	8	4
		white	YES	NO	780013	- 4260348615281	4 x 5	8	4
	6-ply	green	YES	YES	780024	11860852 4260348615298	4 x 5	8	4
				NO	780023	- 4260348615304	4 x 5	8	4
		blue	YES	YES	780002	- 4260348615311	4 x 5	8	4
				NO	780012	- 4260348615328	4 x 5	8	4

Product data sheet

SIZE	PLY	COLOUR	X-RAY	LOOP	REF	PZN (BOX) GTIN	PCS/ BUNDLE	BUNDLE/ BOX	BOXES/ CARTON
20 x 30 cm	4-ply	green	YES	YES	780003	- 4260348615335	4 x 5	5	4
	6-ply	green	YES	YES	780016	- 4260348615342	4 x 5	5	4
				NO	780034	- 4260348615359	4 x 5	5	4
		blue	YES	NO	780036	- 4260348616851	4 x 5	5	4
	white	YES	NO	780037	- 4260348616868	4 x 5	5	4	
35 x 35 cm	8-ply	green	YES	NO	780011	- 4260348615366	4 x 5	10	1
40 x 40 cm	4-ply	green	YES	NO	29107	- 4260348615373	4 x 5	5	4
		blue	YES	YES	780004	- 4260348615380	4 x 5	5	4
		white	YES	YES	780026	- 4260348615397	4 x 5	5	4
	6-ply	green	YES	YES	780005	- 4260348615403	4 x 5	4	4
		blue	YES	YES	780006	- 4260348615410	3 x 6	5	4
				NO	780009	- 4260348615427	3 x 6	5	4
45 x 45 cm	4-ply	green	YES	YES	29175	11859837 4260348615434	4 x 5	5	4
	6-ply	green	YES	YES	250614	11859820 4260348615441	4 x 5	10	1
				NO	250615	- 4260348615458	4 x 5	10	1
		white	YES	YES	930510	11860869 4260348615465	4 x 5	5	4
40 x 80 cm	4-ply	green	YES	NO	29187	- 4260348615472	4 x 5	10	1
50 x 60 cm	4-ply	green	YES	YES	780007	- 4260348615489	4 x 5	5	4
	6-ply	green	YES	YES	780041	- 4260348615496	4 x 5	10	1

Product data sheet

SIZE	PLY	COLOUR	X-RAY	LOOP	REF	PZN (BOX) GTIN	PCS/ BUNDLE	BUNDLE/ BOX	BOXES/ CARTON
NOT pre-washed:									
30 x 45 cm	4-ply	white	YES	NO	780048	- 4260604530280	4 x 5	5	4
28-thread, NOT pre-washed:									
40 x 40 cm	9-ply	blue	YES	NO	780047	- 4260348616677	4 x 5	5	2

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
 The product has to be sterilised using a validated procedure before coming into contact with open wounds.
 The product has not to be used, if the packaging is damaged.
 The product has not to be reused.
 The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: IIa / 7 (annex IX directive 93/42/EEC)
 Conformity assessment: annex VII + V (directive 93/42/EWG)
 The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-09-02

Product description:

20- or 28-thread gauze cloth according to DIN EN 14079,
100 % cotton, with or without loop (30 cm),
pre-washed (unless otherwise stated), with X-RAY contrast chip

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, for haemostasis

Product image (example):



Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE	PLY	PCS/ POUCH	COLOUR	X-RAY	LOOP	REF	PZN (POUCH) GTIN	POUCH/ BOX	BOX/ CTN
8 x 90 cm	4-ply	2	green	YES	NO	530001	- 4260348615533	40	4
					YES	530002	11860415 4260348615540	40	4
			white	YES	YES	540003	- 4260348616615	40	4
	6-ply	2	green	YES	YES	530026	- 4260348616837	28	4
			white	YES	YES	540004	- 4260348615557	28	4
		5	blue	YES	NO	550012	- 4260348615564	20	4
10 x 20 cm	4-ply	10	white	YES	NO	540034	- 4260604532628	15	4
15 x 25 cm	4-ply	2	green	YES	NO	530003	- 4260348615571	40	4

Product data sheet

SIZE	PLY	PCS/ POUCH	COLOUR	X-RAY	LOOP	REF	PZN (POUCH) GTIN	POUCH/ BOX	BOX/ CTN			
20 x 30 cm	4-ply	2	green	YES	YES	530024	- 4260348616608	40	4			
					NO	530004	- 4260348615588	40	4			
		5	white	YES	YES	540005	11860444 4260348615595	40	4			
						530033	- 4260348615601	20	4			
						540006	- 4260348615618	20	4			
		20	white	YES	YES	540033	- 4260348616844	10	4			
	6-ply	5	green	YES	YES	530023	- 4260348615625	20	4			
40 x 40 cm	4-ply	2	green	YES	YES	530006	- 4260348615656	28	4			
						530005	- 4260348615663	12	4			
		5	green	YES	NO	530007	- 4260348615670	10	4			
	YES	530007			- 4260348615670	10	4					
	6-ply	2	green	YES	YES	530014	- 4260348615700	20	4			
						5	blue	YES	NO	550008	- 4260348615724	10
white							YES	YES	540007	- 4260348615731	10	4
45 x 45 cm	4-ply	2	green	YES	YES	530031	- 4260348615748	16	4			
						5	green	YES	YES	530015	11860438 4260348615755	10
		white	YES	YES	540008					11860450 4260348615762	10	4
		2	blue	YES	YES		550013	- 4260348616202	10	4		
		10	blue	YES	NO	550004	- 4260348615779	8	4			
	6-ply	2	green	YES	YES	530008	- 4260348615786	14	4			
						5	green	YES	YES	530009	11860421 4260348615793	9

Product data sheet

SIZE	PLY	PCS/ POUCH	COLOUR	X-RAY	LOOP	REF	PZN (POUCH) GTIN	POUCH/ BOX	BOX/ CTN
45 x 45 cm	6-ply	5	white	YES	YES	540009	11860467 4260348615809	9	4
40 x 80 cm	4-ply	3	green	YES	NO	530013	- 4260348615816	10	4
50 x 60 cm	4-ply	2	white	YES	YES	540010	11860473 4260348615847	12	4
	6-ply	5	green	YES	YES	530021	- 4260348615854	12	4
NOT pre-washed:									
30 x 30 cm	4-ply	10	green	YES	NO	530036	- 4260348615861	10	4
40 x 40 cm	9-ply	5	blue	YES	NO	550005*	- 4260348615878	10	2

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: IIa / 7 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

*Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-09-02

Product description:

cloth made of non-woven fabric (30 g/m²),
without loop,
with X-RAY contrast chip

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, for haemostasis

Product image (example):



Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE	PLY	POUCH	COLOUR	X-RAY	LOOP	REF	PZN (POUCH) GTIN	POUCH/ BOX	BOXES/ CARTON
40 x 60 cm	2-ply	1	white	YES	NO	545001	- 4260604531409	30	4

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: IIa / 7 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-09-02

Product data sheet

Product description:

20-thread bag made of gauze according to DIN EN 14079,
diameter 20 cm or 30 cm, with beaded edge,
100 % cotton,
pre-washed,
available as a flat swab or in pouch form,
with or without loop if required, loop 30 cm,
with X-RAY contrast chip

Field of application / indication:

for short-term storage of the intestine during abdominal surgery (to keep the intestine warm and moist)

Product image (example):



Product variants:

Other variants (size, bag shape, etc.) can be produced according to your requirements.

DIAMETER	PLY	COLOUR	X-RAY	LOOP	REF	GTIN	PCS/BOX
20 cm	6-ply	WHITE	YES	NO	250619	4260348615885	10
23 cm	8-ply	WHITE	YES	YES	250622	4260348615892	10

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: IIa / 7 (annex IX directive 93/42/EEC)

Conformity assessment: annex II without section (4) (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-09-02

Product description:

20-thread bag made of gauze according to DIN EN 14079,
diameter 20 cm or 30 cm,
with beaded edge,
100 % cotton,
pre-washed,
available as a flat swab or in pouch form,
with or without loop if required,
loop 30 cm, with X-RAY contrast chip

Field of application / indication:

for short-term storage of the intestine during abdominal surgery (to keep the intestine warm and moist)

Product image (example):



Product variants:

STERILE sealed (double packed with counting cards for the operating room):

Other variants (size, pouch shape, etc.) can be produced according to your requirements.

DIAMETER	PLY	PCS/ POUCH	COLOUR	X- RAY	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
20 cm without loop	6-ply	1	WHITE	YES	31967	11859866	4260348615915	12	4
23 cm without loop	6-ply	1	WHITE	YES	31973	-	4260348616875	12	4
30 cm (flat) with loop	8-ply	1	WHITE	YES	31972	-	4260348615922	20	4

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-09-02

Product description:

absorbent strip made of 20-thread gauze according to DIN EN 14079,
100 % cotton

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, for haemostasis

Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE		PCS/ POUCH	X-RAY	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
1cm x 30cm	4-ply	1	NO	32467	13568759	4260348616745	100	variable
1cm x 100cm	4-ply	1	NO	32471	11860102	4260348615953	50	variable

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-03-01

Product data sheet

Product description:

17-thread swabs made of gauze according to DIN EN 14079,
 100 % cotton,
 with tucked-in cut edges,
 with a 20 cm long security thread,
 with X-RAY contrast thread

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids, for haemostasis

Product image (example):



Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE	PCS/ POUCH	X-RAY	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
3 x 20 cm	1	YES	233143	-	4260348615984	30	10

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

The security thread is used to locate the compress and not to remove it. The compresses may only be removed from the surgical field using surgical instruments.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-03-90

Product description:

stick made of wood,
head made of wadding (100 % cotton)

diameter:


small head approx. 5 mm

large head approx. 10 mm

Field of application / indication:

for care of secreting wounds, for sucking up blood or other body fluids, only for external use

Product image (example):



Product variants:

SIZE		REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
small head	150 mm of length	2139011	08822428	4260348611252	100	100
	230 mm of length	2139013	08822434	4260348611269	100	100
large head	150 mm of length	2139012	08822440	4260348611276	100	100

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
The product has to be sterilised using a validated procedure before coming into contact with open wounds.
The product has not to be used, if the packaging is damaged.
The product has not to be reused.
The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)
Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)
The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-21-09-05

Product description:

stick made of wood,
head made of wadding (100 % cotton)

diameter:

small head approx. 5 mm

large head approx. 10 mm

Field of application / indication:

for wound treatment, for care of secreting wounds, for sucking up blood or other body fluids,
only for external use

Product image (example):



Product variants:

STERILE sealed (single packaged):

SIZE		REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
150 mm of length	small head	329011	08822457	4260348611283	100 x 1	20
		329012	00050831	4260348611290	100 x 2	20
	large head	329021	08822486	4260348611306	100 x 1	10
		329022	08822492	4260348611313	50 x 2	10

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Is / 4 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCl@ss:

34-21-09-05

Product description:

white fabric made of 45 % viscose and 55 % polyester,
elongation approx. 50 %,
light, soft,
air-permeable,
selvedges

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, only for external use

Product image (example):



Product variants:

naked, box à 20, 50 or 100 bandages

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOX/ CARTON
4cm x 4m	1050501	08821914	4260348611320	20	48
6cm x 4m	1050502	08821920	4260348611337	20	32

Product data sheet

8cm x 4m	1050503	08821937	4260348611344	20	24
10cm x 4m	1050504	08821943	4260348611351	20	20
12cm x 4m	1050505	08821966	4260348611368	20	16
4cm x 4m	1050516	10180090	4260348611375	50	24
6cm x 4m	1050517	10180109	4260348611382	50	16
8cm x 4m	1050518	10180115	4260348611399	50	12
10cm x 4m	1050519	10180121	4260348611405	50	8
4cm x 4m	1050510	04793954	4260348611412	100	12
6cm x 4m	1050511	08821972	4260348611429	100	8
8cm x 4m	1050512	08821989	4260348611436	100	6
10cm x 4m	1050513	08821995	4260348611443	100	4

individually cellophane wrapped:



SIZE	REF	PZN (POUCH WITH 1 PIECE)	GTIN	PCS/ POUCH	POUCH/ BOX	BOX/ CTN
4cm x 4m	25102304	12502735	4260348617872	1	20	48
6cm x 4m	25102306	12502741	4260348617889	1	20	32
8cm x 4m	25102308	12502758	4260348617896	1	20	24
10cm x 4m	25102310	12502764	4260348617902	1	20	20

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-02-01

Product description:

white fabric made of 45 % viscose and 55 % polyester,
 elongation approx. 50 %,
 light, soft,
 air-permeable,
 selvedges

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, only for external use

Product variants:

STERILE sealed (single packaged):

SIZE	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
4cm x 4m	32517	11860183	4260348611498	50 x 1	8
6cm x 4m	32518	11860208	4260348611504	50 x 1	8
8cm x 4m	32525	11860220	4260348611511	50 x 1	8
10cm x 4m	32532	11860289	4260348611528	50 x 1	8
12cm x 4m	32520	11860214	4260348611535	30 x 1	8

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
 The product has not to be used, if the packaging is damaged.
 The product has not to be reused.
 The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-02-01

Product description:

white fabric made of 45 % viscose and 55 % polyester,
 elongation approx. 50 %,
 light, soft,
 air-permeable,
 selvedges

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, only for external use

Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE	POUCH	REF	PZN (POUCH)	GTIN	POUCH/BOX	BOXES/CARTON
4cm x 4m	1	32513	-	4260348616011	30	6
6cm x 4m	1	32475	-	4260348616028	20	7
8cm x 4m	1	32476	-	4260348616035	60	4
10cm x 4m	1	32454	-	4260348616042	50	4
12cm x 4m	1	32453	-	4260348616059	40	4

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
 The product has not to be used, if the packaging is damaged.
 The product has not to be reused.
 The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-02-01

Product data sheet

Product description:

white fabric made of 70 % viscose and 30 % polyamide,
 elongation approx. 100 %,
 with crimped polyamide fibres,
 light, soft,
 air-permeable,
 selvedges

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, only for external use

Product image (example):



Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
4cm x 4m	1050601	04793931	4260348611542	20	24
6cm x 4m	1050602	03842900	4260348611559	20	16
8cm x 4m	1050603	03842946	4260348611566	20	12
10cm x 4m	1050604	03842981	4260348611573	20	10
12cm x 4m	1050605	04793948	4260348611580	20	8
4cm x 4m	1050611	04793894	4260348611597	100	8
6cm x 4m	1050612	04793902	4260348611603	100	8
8cm x 4m	1050613	04793919	4260348611610	100	6
10cm x 4m	1050614	04793925	4260348611627	100	4

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-02-01

Product description:

white fabric made of 20-thread gauze according to DIN EN 14079,
100 % cotton,
no stretch,
light, skin-friendly, soft,
air-permeable,
selvedges

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, only for external use

Product image (example):



Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
4cm x 4m	1739011	04970254	4260348611634	20	60
6cm x 4m	1739013	04868161	4260348611641	20	60
8cm x 4m	1739015	04868178	4260348611658	20	48
10cm x 4m	1739017	04868215	4260348611665	20	36
12cm x 4m	1739019	04868221	4260348611672	20	36

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-01-01

Product description:

raw white fabric made of 54 % cotton, 37 % polyester,
 8,5 % spandex and 0,5 % polyamide,
 elongation approx. 100 %,
 medium to short stretch,
 permanently elastic,
 selvedges

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints,
 only for external use

Product image (example):

Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
6cm x 5m	815067	07200274	4260348611689	10	40
8cm x 5m	815087	07200268	4260348611696	10	28
10cm x 5m	815107	07200251	4260348611702	10	24
12cm x 5m	815127	07200245	4260348611719	10	20
15cm x 5m	815158	07200096	4260348611726	10	16
20cm x 5m	815207	07200073	4260348611733	10	12

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
The product has to be sterilised using a validated procedure before coming into contact with open wounds.
The product has not to be used, if the packaging is damaged.
The product has not to be reused.
The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)
Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)
The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-04-01

Product description:

raw white fabric made of 54 % cotton, 37 % polyester, 8,5 % spandex and 0,5 % polyamide,
elongation approx. 100 %, medium to short stretch,
permanently elastic, selvages

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints,
only for external use

Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE	POUCH	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
12cm x 5m	1	32460	-	4260604533731	15	8

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
The product has not to be used, if the packaging is damaged.
The product has not to be reused.
The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-04-01

Product description:

white fabric made of 60 % cotton, 37 % polyester, 1 % polyamide and 2 % elastane, elongation approx. 90 %, permanently elastic, selvages

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints, for mild compression, only for external use

Product image (example):

Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
6 cm	1725815001	08822061	4260348611740	10	40
8 cm	1725815002	08822078	4260348611757	10	28
10 cm	1725815003	08822084	4260348611764	10	24
12 cm	1725815004	08822090	4260348611771	10	20
15 cm	1725815006	11859671	4260348611788	10	16
20 cm	1725815007	11859688	4260348611795	10	12

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-04-01

Product description:

white fabric made of 60 % cotton, 39 % polyamide and 1 % elastane,
 elongation approx. 80 %,
 permanently elastic,
 selvages,
 with bandage clips

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints,
 for mild compression, only for external use

Product image (example):

Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
6cm x 5m	815066	04794741	4260348611801	10	10
8cm x 5m	815086	04794758	4260348611818	10	10
10cm x 5m	815106	04794764	4260348611825	10	10
12cm x 5m	815126	04794770	4260348611832	10	5

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: IX (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-03-01

Product description:

white fabric made of 60 % cotton, 39 % polyamide and 1 % elastane,
 elongation approx. 92 %,
 permanently elastic,
 selvedges,
 with bandage clips

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints,
 for mild compression, only for external use

Product variants:

STERILE sealed (single packaged):

SIZE	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
6cm x 5m	32498	11860160	4260348611849	18 x 1	8
8cm x 5m	32491	11860125	4260348611856	17 x 1	8
10cm x 5m	32482	11860119	4260348611863	20 x 1	6
12cm x 5m	32457	11860065	4260348611870	20 x 1	8
15cm x 5m	32499	11860177	4260348611887	22 x 1	3

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
 The product has not to be used, if the packaging is damaged.
 The product has not to be reused.
 The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCl@ss:

34-19-03-01

Product description:

white fabric made of 60 % cotton, 39 % polyamide and 1 % elastane,
elongation approx. 92 %,
permanently elastic,
selvedges,
with bandage clips

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints,
for mild compression, only for external use

Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE	POUCH	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
6cm x 5m	1	32497	-	4260348616066	20	8
8cm x 5m	1	32470	-	4260348616073	15	8
10cm x 5m	1	32489	-	4260348616080	15	8
12cm x 5m	1	32490	-	4260348616097	15	8
15cm x 5m	1	32529	-	4260348616103	13	8

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
The product has not to be used, if the packaging is damaged.
The product has not to be reused.
The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCl@ss:

34-19-03-01

Product description:

fabric made of 60 % cotton, 39 % polyamide and 1 % elastane,
 elongation approx. 90 %,

permanently elastic,

selvedges,

with bandage clips

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints,
 for mild compression, only for external use

Product image (example):

Product variants:

SIZE	COLOUR	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
6cm x 5m	RED	815061	04794787	4260348611894	10	10
	BLUE	815062	04794793	4260348611900	10	10
8cm x 5m	RED	815081	04794824	4260348611917	10	10
	BLUE	815082	04794830	4260348611924	10	10
10cm x 5m	RED	815101	04794913	4260348611931	10	10
	BLUE	815102	04794936	4260348611948	10	10

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
The product has to be sterilised using a validated procedure before coming into contact with open wounds.
The product has not to be used, if the packaging is damaged.
The product has not to be reused.
The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)
Conformity assessment: IX (regulation (EU) 2017/745)
The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-03-01

Product data sheet

Product description:

fabric made of 88 % cotton, 11 % polyamide and 1 % elastane,
stretch approx. 80 %, high cotton content,
skin-friendly, breathable,
permanently elastic, selvedges,
with bandage clips

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints,
for mild compression, only for external use

Product image (example):



Product variants:

SIZE	COLOUR	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
6cm x 5m	RED	11440650	03759470	4260348611955	10	10
	BLUE	11450650	03759487	4260348611962	10	10
	YELLOW	11470650	03759599	4260348611979	10	10
	PURPLE	11480650	03759671	4260348611986	10	10

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: IX (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-03-01

Product description:

white fabric made of 100 % cotton,
elongation approx. 90 %,
textile elastic,
loop-edged,
according to DIN 61632

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints,
for mild compression, only for external use

Product image (example):



Product variants:

naked (without bandage clips):

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
4cm x 5m	813042	04881345	4260348611993	10	30
6cm x 5m	813062	08822115	4260348612006	10	20
8cm x 5m	813082	08822121	4260348612013	10	15
10cm x 5m	813102	08822138	4260348612020	10	15
12cm x 5m	813122	08822144	4260348612037	10	15
15cm x 5m	813151	08822150	4260348612044	10	6
20cm x 5m	813202	02898235	4260348612051	10	8
30cm x 5m	813302	07200541	4260348612075	10	16

individually cellophaned (with bandage clips):

SIZE	REF	PZN (PCS)	GTIN	PCS/BOX	BOXES/CARTON
6cm x 5m	813063	08822173	4260348612082	10	10
8cm x 5m	813083	08822196	4260348612099	10	10
10cm x 5m	813103	08822204	4260348612105	10	5
12cm x 5m	813123	08822210	4260348612112	10	5
30cm x 5m	813303	07200558	4260348612068	1	variable

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-03-01

Product description:

white fabric made of 100 % cotton,
 elongation approx. 90 %,
 textile elastic,
 loop-edged,
 according to DIN 61632

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints,
 for mild compression, only for external use

Product variants:

STERILE sealed (single packaged):

SIZE	POUCH	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
6cm x 5m	1	32527	11860237	4260348612129	18	8
8cm x 5m	1	32494	11860131	4260348612136	18	8
10cm x 5m	1	32495	11860148	4260348612143	20	6
12cm x 5m	1	32496	11860154	4260348612150	20	6
15cm x 5m	1	32528	11860243	4260348612167	12	6

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
 The product has not to be used, if the packaging is damaged.
 The product has not to be reused.
 The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCl@ss:

34-19-03-01

Product description:

white fabric made of 100 % cotton,
 elongation approx. 90 %,
 textile elastic,
 loop-edged,
 according to DIN 61632

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints,
 for mild compression, only for external use

Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE	POUCH	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
6cm x 5m	1	32533	-	4260348616110	15	8
8cm x 5m	1	32493	-	4260348616127	15	8
10cm x 5m	1	32492	-	4260348616134	30	3
12cm x 5m	1	32469	-	4260348616141	30	3
15cm x 5m	1	32505	-	4260348616158	10	6

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
 The product has not to be used, if the packaging is damaged.
 The product has not to be reused.
 The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCl@ss:

34-19-03-01

Product description:

skin-coloured fabric made of 81 % cotton, 16 % polyamide and 3 % elastane,
elongation approx. 165 %, permanently elastic,
with bandage clips

Field of application / indication:

for compression, for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints, only for external use

Product image (example):



Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
6cm x 7m	210266	04793003	4260348612242	8	10
8cm x 7m	210267	04793026	4260348612181	8	10
10cm x 7m	210268	04793032	4260348612198	8	5
12cm x 7m	210269	04793049	4260348612204	8	5

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: IX (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-05-01

Product description:

skin-coloured fabric made of 100 % cotton,
 elongation approx. 80 %,
 textile elastic,
 with bandage clips

Field of application / indication:

for compression, for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints, only for external use

Product image (example):

Product variants:

SIZE	REF	PZN (BOX)	GTIN	POUCH/ BOX	BOXES/ CARTON
6cm x 5m	24126206509	04793055	4260348612259	10	16
8cm x 5m	24126208509	04793061	4260348612266	10	12
	24126208510	09082011	4260348612273	80	-
10cm x 5m	24126210509	04793078	4260348612280	10	10
	24126210510	09082005	4260348612297	60	-
12cm x 5m	24126212509	04793084	4260348612303	10	10
	24126212510	09081997	4260348612310	60	-

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-03-01

Product description:

skin-coloured fabric made of 100 % cotton, elongation approx. 80 %, textile elastic, with bandage clips

Field of application / indication:

for compression, for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints, only for external use

Product variants:

STERILE sealed (double packed with counting cards for the operating room):

SIZE	POUCH	REF	PZN (POUCH)	GTIN	POUCH/ BOX	BOXES/ CARTON
6cm x 5m	1	32530	11860272	4260348616165	12	8
8cm x 5m	1	32464	11860088	4260348616172	12	8
10cm x 5m	1	32462	11860071	4260348616189	15	6
12cm x 5m	1	32455	11860059	4260348616196	12	6

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-03-01

Product description:

compression bandage,
consisting of 2 short-stretch bandages for the classic opposing bandage technique,
skin-coloured fabric made of 100 % cotton,
elongation approx. 80 %,
textile elastic,
with bandage clips

Field of application / indication:

for compression, for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints, only for external use

Product image (example):



Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
bandage 1: 8cm x 5m bandage 2: 10cm x 5m	24140	12396954	4260348614000	2	64

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-09-01

Product description:

skin-coloured fabric made of 67 % cotton and 33 % polyamide,
 elongation approx. 40 %,

permanently elastic,

particularly strong compression,

selvages,

with bandage clips

Field of application / indication:

for compression, for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints, only for external use

Product image (example):

Product variants:

SIZE	REF	PZN	GTIN	PCS/BOX	BOXES/CARTON
6cm x 5m	210231	04793227	4260348612327	10	10
8cm x 5m	210232	04793233	4260348612334	10	10
10cm x 5m	210233	04793262	4260348612341	10	5
12cm x 5m	210234	04793339	4260348612358	10	5

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: IX (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-03-01

Product description:

bonded non-woven fibre product with elastane threads,
 consisting of approx. 50 % non-woven fabric (polypropylene),
 approx. 10 % elastane and approx. 40 % synthetic adhesive,
 elongation approx. 100 %, only adheres to itself,
 removable without residue, elastic,
 can be torn by hand, water/oil-resistant,
 breathable, individually cellophane-wrapped

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints,
 only for external use

Product image (example):



Product variants:

SIZE	COLOUR	REF	PZN (PCS)	GTIN	PCS/BAG	BAGS/BOX	BOXES/CARTON
2,5cm x 4,5m	blue	107004	12144750	4260348616387	1	20	18
5cm x 4,5m	blue	107001	12535893	4260604532598	1	10	30
	red	107007	12535918	4260348618183	1	10	30

Product data sheet

7,5cm x 4,5m	blue	107006	12535901	4260348618176	1	10	24
	red	107008	12535924	4260348618190	1	10	24
10cm x 4,5m	blue	107003	15436401	4260604532567	1	10	18
	red	107009	-	4260604533793	1	10	18

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-02-02

Product description:

white fabric, elongation approx. 100 %, only adheres to itself

LATEX CONTENT: 50 % cotton and 50 % polyamide, latex

LATEX-FREE: 46 % cotton and 54 % polyamide

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints, only for external use

Product image (example):



Product variants:

	SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
WITH LATEX	4cm x 20m	1061101	11152635	4260348612365	1	20
	6cm x 20m	1061102	04807136	4260348612372	1	20
	8cm x 20m	1061103	04807142	4260348612389	1	16
	10cm x 20m	1061104	04807159	4260348612396	1	12
	12cm x 20m	1061105	04807165	4260348612402	1	12
LATEXFREE	4cm x 20m	1061401	12397830	4260348617766	1	24
	6cm x 20m	1061402	12397847	4260348617773	1	18
	8cm x 20m	1061403	12397876	4260348617780	1	12
	10cm x 20m	1061404	12397882	4260348617797	1	12
	12cm x 20m	1061405	12397899	4260348617803	1	10

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-02-02

Product description:

fabric made of 39 % cotton and 61 % polyamide,
elongation approx. 60 %,
with a latex coating,
therefore only adheres to itself

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints,
only for external use

Product image (example):



Product variants:

SIZE	COLOUR	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
4cm x 20m	RED	10511011	04793977	4260348612419	1	36
	BLUE	10511012	04793960	4260348612426	1	36
6cm x 20m	RED	10511021	03842797	4260348612433	1	36
	BLUE	10511022	03842840	4260348612440	1	36

Product data sheet

8cm x 20m	RED	10511031	03759725	4260348612457	1	36
	BLUE	10511032	03759731	4260348612464	1	36
	GELB	10511033	03759748	4260348612471	1	36
	GREEN	10511034	03759837	4260348612488	1	36
10cm x 20m	RED	10511041	03759984	4260348612495	1	18
	BLUE	10511042	03759990	4260348612501	1	18

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: IX (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-02-02

Product description:

white fabric, elongation approx. 100 %, only adheres to itself

LATEX CONTENT: 42 % cotton, 29 % polyamide and 29 % rayon, latex

LATEX-FREE: 46 % cotton, 27 % polyamide and 27 % rayon

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints, only for external use

Product image (example):



Product variants:

4 metres long:

	SIZE	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
WITH LATEX	4cm x 4m	1061201	04794652	4260348612518	1	84
	6cm x 4m	1061202	04794669	4260348612525	1	54
	8cm x 4m	1061203	04794675	4260348612532	1	48
	10cm x 4m	1061204	04794681	4260348612549	1	36
	12cm x 4m	1061205	04794698	4260348612556	1	26
LATEXFREE	4cm x 4m	1061301	12397787	4260348617285	1	84
	6cm x 4m	1061302	12397793	4260348617308	1	54
	8cm x 4m	1061303	12397801	4260348617315	1	48
	10cm x 4m	1061304	12397818	4260348617322	1	36
	12cm x 4m	1061305	12397824	4260348617612	1	26

20 metres long:

	SIZE	REF	PZN (BOX)	GTIN	PCS/BOX	BOX/CTN
WITH LATEX	4cm x 20m	1051101	08822227	4260348612563	1	20
	6cm x 20m	1051102	08822233	4260348612570	1	20
	8cm x 20m	1051103	08822256	4260348612587	1	16
	10cm x 20m	1051104	08822262	4260348612594	1	12
	12cm x 20m	1051105	08822279	4260348612600	1	12
LATEXFREE	4cm x 20m	1051201	12397735	4260348617711	1	22
	6cm x 20m	1051202	12397741	4260348617728	1	16
	8cm x 20m	1051203	12397758	4260348617735	1	12
	10cm x 20m	1051204	12397764	4260348617742	1	10
	12cm x 20m	1051205	12397770	4260348617759	1	8

length specifications in a stretched state

Warnings / contraindications:

- The product is not to be used in the case of an existing allergy to the material.
- The product has to be sterilised using a validated procedure before coming into contact with open wounds.
- The product has not to be used, if the packaging is damaged.
- The product has not to be reused.
- The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

- Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)
- Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)
- The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-02-02

Product description:

fabric made of 72 % cotton and 28 % polyamide,
elongation approx. 60 %,
with a latex coating,
therefore only adheres to itself

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints,
only for external use

Product image (example):



Product variants:

SIZE	COLOUR	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
4cm x 20m	RED	10611011	04794190	4260348612617	1	22
	BLUE	10611012	04794238	4260348612624	1	22
	YELLOW	10611013	11859429	4260348612631	1	22
	GREEN	10611014	11859435	4260348612648	1	22
6cm x 20m	RED	10611021	04794422	4260348612655	1	16
	BLUE	10611022	04808035	4260348612662	1	16
	YELLOW	10611023	11859441	4260348612679	1	16
	GREEN	10611024	11859458	4260348612686	1	16

SIZE	COLOUR	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
8cm x 20m	RED	10611031	04794468	4260348612693	1	12
	BLUE	10611032	04794474	4260348612709	1	12
	YELLOW	10611033	11859464	4260348612716	1	12
	GREEN	10611034	11859470	4260348612723	1	12
10cm x 20m	RED	10611041	04794592	4260348612730	1	10
	BLUE	10611042	04794600	4260348612747	1	10

length specifications in a stretched state

Warnings / contraindications:

- The product is not to be used in the case of an existing allergy to the material.
- The product has to be sterilised using a validated procedure before coming into contact with open wounds.
- The product has not to be used, if the packaging is damaged.
- The product has not to be reused.
- The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

- Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)
- Conformity assessment: IX (regulation (EU) 2017/745)
- The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-02-02

Product description:

white fabric made of 60 % cotton, 39 % polyamide and 1 % elastane,
elongation approx. 80 %,
with a latex coating,
therefore only adheres to itself

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints,
for mild compression, only for external use

Product image (example):



Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
6cm x 5m	820060	04794706	4260348612754	1	64
8cm x 5m	820080	04794712	4260348612761	1	64
10cm x 5m	820100	04794729	4260348612778	1	32
12cm x 5m	820120	04794735	4260348612785	1	32

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: IX (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-03-02

Product data sheet

Product description:

skin-coloured fabric made of 60 % cotton, 39 % polyamide and 1 % elastane,
 elongation approx. 90 %,
 with a latex coating,
 therefore only adheres to itself

Field of application / indication:

for compression, for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints, only for external use

Product image (example):



Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
6cm x 5m	2481120652	00537036	4260348612792	1	64
8cm x 5m	2481120852	00537042	4260348612808	1	64
10cm x 5m	2481121052	00537059	4260348612808	1	32
12cm x 5m	24811221252	00537065	4260348612822	1	32

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: IX (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-03-02

Product description:

skin-coloured fabric made of 100 % cotton,
elongation approx. 60 %,
with particularly skin-friendly acrylate adhesive,
textile-elastic,
selvedges

Field of application / indication:

for compression, to support and relieve joints, only for external use

Product image (example):



Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
6cm x 2,5m	100003	04793345	4260348612839	1	72
8cm x 2,5m	100002	04793405	4260348612846	1	72
10cm x 2,5m	100001	04793411	4260348612853	1	36

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: IX (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-06-01

Product description:

skin-coloured fabric made of 100 % cotton,
bi-elastic approx. 30 % across,
with particularly skin-friendly acrylate adhesive,
textile-elastic,
selvedges

Field of application / indication:

for compression, for the fixation of wound dressings, splints, cannulas or similar medical devices, to support and relieve joints, only for external use

Product image (example):



Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
6cm x 2,5m	210331	04793428	4260348612860	1	72
8cm x 2,5m	210332	04793440	4260348612877	1	72
10cm x 2,5m	210333	04793463	4260348612884	1	36

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: IX (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-06-01

Product description:

100 % polyester,
85 g/m²,
white,
breathable,
water-repellent properties,
tearable

Field of application / indication:

for padding
(plaster-, cast- and compression-bandages), only for external use

Product image (example):



Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
4 cm x 3 m	420025	11058681	4260348612891	12	40
6 cm x 3 m	420026	11058698	4260348612907	12	28
7,5 cm x 2,7 m	420002	13831464	4260348616776	6	-
10 cm x 3 m	420035	15881153	4260604532727	6	30
10 cm x 3 m	420027	11058706	4260348612914	12	16
15 cm x 3 m	420028	11058712	4260348612921	12	12

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-20-03-01

Product data sheet

Product description:

raw white fabric made of 98 % cotton and 2 % elastane,
 elongation approx. 130 %,
 permanently elastic,
 breathable,
 skin-friendly,
 supple and soft during use

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices (if a padding or a soft and gentle application is important), only for external use

Product image (example):



Product variants:

STERILE sealed (single packaged):

SIZE	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
10cm x 6m	317011	10763609	4260348612938	20 x 1	6

length specifications in a stretched state

Warnings / contraindications:

- The product is not to be used in the case of an existing allergy to the material.
- The product has not to be used, if the packaging is damaged.
- The product has not to be reused.
- The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Sterile procedure pack according to Art. 12 (3) of Directive 93/42/EEC

Conformity assessment: -

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-02-01

Product description:

white fabric made of 100 % cotton,
supple and soft,
great cushioning effect,
high absorbency

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices
(if a padding or a soft and gentle application is important), only for external use

Product image (example):



Product variants:

SIZE	REF	PZN	GTIN	PCS/ BAG	BAGS/ CARTON
11,4cm x 3,7m	317015	-	4260604532673	10	10
11,4cm x 3,7m	317017	-	4260604532697	100	1

length specifications in a stretched state

Warnings / contraindications:

- The product is not to be used in the case of an existing allergy to the material.
- The product has to be sterilised using a validated procedure before coming into contact with open wounds.
- The product has not to be used, if the packaging is damaged.
- The product has not to be reused.
- The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-02-01

Product description:

white fabric made of 100 % cotton,
supple and soft,
great cushioning effect,
high absorbency

Field of application / indication:

for the fixation of wound dressings, splints, cannulas or similar medical devices
(if a padding or a soft and gentle application is important), only for external use

Product variants:

STERILE sealed (single packaged):

SIZE	POUCH	REF	PZN	GTIN	POUCH/ BOX	BOXES/ CARTON
11,4cm x 3,7m	1	317016	-	4260604532758	50	4

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
The product has not to be used, if the packaging is damaged.
The product has not to be reused.
The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Is / 4 (annex IX directive 93/42/EEC)
Conformity assessment: annex VII + V (directive 93/42/EWG)
The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-02-01

Product description:

fabric made of rayon and polyamide,
 elongation approx. 50 to 60 %,
 low resting pressure,
 high working pressure,
 high layer adhesion,
 stable edges,
 skin-friendly,
 easy to put on and take off

Field of application / indication:

Zinc adhesive bandages are used as compression bandages for phlebological conditions such as chronic venous insufficiency including venous leg ulcers, venous oedema, thrombophlebitis and phlebothrombosis. Semi-rigid bandages are used as support and relief bandages for injuries to the musculoskeletal system, e.g. sprains, dislocations, strains, bruises, tenosynovitis, further treatment after plaster casts, to reduce post-traumatic oedema., only for external use

Product image (example):



Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
8cm x 7m	490002	12528195	4260348618213	10	10
10cm x 5m	490003	12528203	4260348618220	10	10
10cm x 7m	490004	12528232	4260348618237	10	10
10cm x 10m	490005	12528249	4260348618244	10	10

length specifications in a stretched state

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product should not be used in peripheral arterial occlusive diseases, decompensated cardiac insufficiency, septic phlebitis, phlegmasia coerulea dolens.

The product is not to be used in contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be stored horizontally, protected from sunlight, dry and at a maximum temperature of 25°C.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-20-01-01

Product description:

aluminium bandage clips with white rubber cord

Field of application / indication:

for fixation of bandages, only for external use

Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
approx. 1,5 x 3 cm	2641010	11224323	4260348612945	100	64

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: IX (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-10-03

Product description:

white fabric made of 73 % polyamide and 27 % polyurethane,
 circular weave, wide-meshed,
 stretchy, pressure- and crease-free application,
 stitch-resistant, can be cut at any point without tearing the cuts

Field of application / indication:

for the fixation of wound dressings, only for external use

Product image (example):

Product variants:

SIZE (UNSTRETCHED)	FIELD OF APPLICATION	REF	PZN (BOX)	GTIN	ROLLS/ BOX	BOXES/ CARTON
SIZE 1 1,4cm x 25m	finger	1050822	08822411	4260348612952	1	57
SIZE 2 1,8cm x 25m	hand, arm, foot	1050824	08822405	4260348612969	1	57
SIZE 3 2,3cm x 25m	leg	1050825	08822397	4260348612976	1	57
SIZE 4 3,5cm x 25m	head	1050827	08822380	4260348612983	1	57
SIZE 5 5,0cm x 25m	small torso, shoulder	1050829	08822374	4260348612990	1	30
SIZE 6 5,8cm x 25m	large torso, shoulder	1050830	08822368	4260348613003	1	30

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-08-03

Product data sheet

Product description:

white fabric made of 100 % cotton,
 bleached without chlorine,
 circular knitted gauze bandage, which is very elastic due to the knitting method,
 can be stretched to three times its width,
 becomes tighter when stretched,
 easy to put on,
 wrinkle-free fit,
 non-slip,
 air-permeable

Field of application / indication:

for the fixation of wound dressings, only for external use

Product image (example):



Product variants:

SIZE (UNSTRETCHED)	FIELD OF APPLICATION	REF	PZN (BOX)	GTIN	ROLLS/ BOX	BOXES/ CARTON
SIZE 1 1,5cm x 20m	finger and toe, small	1050731	08822285	4260348613010	1	64
SIZE 2 2,5cm x 20m	finger and toe, large	1050732	08822291	4260348613027	1	33

Product data sheet

SIZE (UNSTRETCHED)	FIELD OF APPLICATION	REF	PZN (BOX)	GTIN	ROLLS/ BOX	BOXES/ CARTON
SIZE 4 4cm x 20m	hand and arm	1050733	08822316	4260348613034	1	29
SIZE 6 6cm x 20m	arm, leg, head	1050734	08822322	4260348613041	1	60
SIZE 8 8cm x 20m	head, leg, thigh	1050735	08822339	4260348613058	1	48
SIZE 10 10cm x 20m	head, small torso	1050736	08822345	4260348613065	1	39
SIZE 12 12cm x 20m	head, medium sized torso	1050738	14398294	4260604532468	1	24
SIZE 16 16cm x 20m	large torso	1050737	8822351	4260348613072	1	18

Warnings / contraindications:

- The product is not to be used in the case of an existing allergy to the material.
- The product has to be sterilised using a validated procedure before coming into contact with open wounds.
- The product has not to be used, if the packaging is damaged.
- The product has not to be reused.
- The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

- Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)
- Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)
- The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-08-04

Product data sheet

Product description:

white fabric made of 67 % cotton and 33 % rayon,
 seamless knitted tubular bandage with high longitudinal and transverse elasticity,
 measured, tied and rolled ready for use,
 quick and easy to apply,
 wrinkle-free fit,
 non-slip,
 air-permeable

Field of application / indication:

for the fixation of wound dressings, only for external use

Product image (example):



Product variants:

DIAMETER (UNSTRECHTED)	APPLICATION	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
Size 1: 2,5 cm	finger (adults)	1050739	00007189	4260348613089	40	30

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-08-04


Product description:

white fabric made of 59 % viscose,
41 % dorlastan,
circular knitted bandage,
lengthwise and crosswise elastic,
easy to put on, wrinkle-free fit,
non-slip, air-permeable,
woven-in colour stripe as size indicator

Field of application / indication:

for the fixation of wound dressings, only for external use

Product image (example):



Product variants:

SIZE (UNSTRETCHED)	FIELD OF APPLICATION (PERIMETER)	REF	PZN (BOX)	GTIN	ROLLS/ BOX	BOXES/ CARTON
SIZE 1 3,5cm x 10m	child's arm, foot (9 – 18 cm)	1060800	12407824	4260348616493	1	35
SIZE 2 5,5cm x 10m	arm, child's leg (14 – 24 cm)	1060801	12407830	4260348616509	1	28
SIZE 3 8cm x 10m	child's head, arm (24 – 40 cm)	1060802	12407847	4260348616516	1	20
SIZE 4 12cm x 10m	head, torso, leg (35 – 64 cm)	1060803	12407853	4260348616523	1	14
SIZE 5 20cm x 10m	torso (64 – 130 cm)	1060804	12407876	4260348616530	1	7

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-19-08-04

Product description:

non-elastic tape bandage made of 100 % cotton,
skin-friendly hotmelt adhesive,
easy to tear thanks to the serrated edge

Field of application / indication:

for stabilizing joints, muscles, and ligaments, to protect them against overstrain and to reduce the risk of injury and to support sprains, only for external use

Product image (example):



Product variants:

SIZE	COLOUR	REF	PZN (BOX)	GTIN	ROLLS/ BOX	BOXES/ CARTON
2 cm x 10 m	white	1672101	01561108	4260348613096	12	10
		1672111	04881813	4260348613102	1	variable
2,5 cm x 10 m	white	1672109	11859524	4260348613119	12	10
3,8 cm x 10 m	white	1672102	01566181	4260348613126	12	10
		1672112	04881486	4260348613133	1	variable
3,8 cm x 10 m	red	1672103	01566755	4260348613140	12	10
		1672113	04881457	4260348613157	1	variable
3,8 cm x 10 m	blue	1672104	01566904	4260348613164	12	10
		1672114	04881434	4260348613171	1	variable

Product data sheet

SIZE	COLOUR	REF	PZN (BOX)	GTIN	ROLLS/ BOX	BOXES/ CARTON
3,8 cm x 10 m	yellow	1672105	01574140	4260348613188	12	10
		1672115	04881397	4260348613195	1	variable
3,8 cm x 10 m	green	1672106	04881492	4260348613201	12	10
		1672116	04880050	4260348613218	1	variable
3,8 cm x 10 m	black	1672118	15241229	4260604532536	12	10
		1672119	15241212	4260604532543	1	variable
5 cm x 10 m	white	1672107	04879070	4260348613225	6	10
		1672117	04878030	4260348613232	1	variable

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:


34-20-01-02

Product description:

elastic special plaster made of 90 % cotton and 10 % elastic urethane fibres,
80 % stretchable (lengthwise),
skin-friendly acrylate adhesive applied in waves,
permeable to air and liquids but water-resistant and quick-drying,
comfortable to wear even over several days,
back with grid for easy cutting to size

Field of application / indication:

for the protection and stabilization of the muscle and joint functions during athletic stress and in the context of orthopedic treatment, for kinesiological therapy, only for external use

Product image (example):

Product variants:

SIZE	COLOUR	REF	PZN (BOX)	GTIN	ROLLS/ BOX	BOXES/ CARTON
5 cm x 5 m (unstretched)	red	1672200	09303653	4260348613249	1	variable
		1672210	09303676	4260348613256	12	12
	blue	1672201	09303682	4260348613263	1	variable
		1672211	09303699	4260348613270	12	12

SIZE	COLOUR	REF	PZN (BOX)	GTIN	ROLLS/ BOX	BOXES/ CARTON
5 cm x 5 m (unstretched)	skin-colour	1672202	09303707	4260348613287	1	variable
		1672212	09303713	4260348613294	12	12
	black	1672203	09303736	4260348613300	1	variable
		1672213	09303742	4260348613317	12	12
	white	1672208	10830228	4260348613324	1	variable
		1672218	10830234	4260348613331	12	12
	orange	1672207	10830205	4260348613348	1	variable
		1672217	10830211	4260348613355	12	12
	mint	1672206	10830180	4260348613362	1	variable
		1672216	10830197	4260348613379	12	12
	yellow	1672223	11286146	4260348613386	1	variable
		1672219	11286123	4260348613393	12	12

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-20-01-02

Product description:

white non-woven (viscose, polyester) with skin-friendly hotmelt adhesive,
permeable to air and water vapour,
with non-adhesive central,
highly absorbent wound-pad (PE, viscose and polyester)

Field of application / indication:

for covering and protecting the wound from infection, primary or secondary dressing over acute wounds and chronic wounds, for post-operative wound care, only for external use

Product image (example):



Product variants:

STERILE sealed (single packaged):

SIZE	SIZE WOUND PAD	REF	PZN (BOX)	GTIN	PCS/BOXES	BOXES/CARTON
5 cm x 7,2 cm	2,5 cm x 4 cm	718000	00466827	4260348613409	100 x 1	24
6 cm x 10 cm	3 cm x 6 cm	718013	10821531	4260348613423	50 x 1	24
8 cm x 10 cm	4 cm x 5 cm	718003	00466862	4260348613430	50 x 1	24
8 cm x 15 cm	4 cm x 10 cm	718005	00466891	4260348613447	50 x 1	24
10 cm x 20 cm	5 cm x 15 cm	718006	00466916	4260348613454	50 x 1	24
10 cm x 25 cm	5 cm x 20 cm	718007	10821548	4260348613461	50 x 1	12
10 cm x 30 cm	5 cm x 25 cm	718008	10821554	4260348613478	50 x 1	12

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Is / 4 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-18-03-03

Product description:

transparent, semi-permeable PU film (polyurethane) with skin-friendly acrylate adhesive,
permeable to air and water vapour,
waterproof and water-tight,
with non-adhesive central, highly absorbent wound-pad

Field of application / indication:

for covering and protecting the wound from infection, primary or secondary dressing over acute wounds and chronic wounds, for post-operative wound care, only for external use

Product image (example):



Product variants:

STERILE sealed (single packaged):

SIZE	SIZE WOUND PAD	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
5 cm x 7 cm	2,5 cm x 4 cm	719005	10821560	4260348613485	50 x 1	24
8 cm x 10 cm	4 cm x 5 cm	719007	10821577	4260348613492	50 x 1	24
10 cm x 15 cm	5 cm x 10 cm	719009	10821637	4260348613508	50 x 1	24
10 cm x 30 cm	5 cm x 25 cm	719011	10821643	4260348613515	50 x 1	12

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

sterilised

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: Is / 4 (annex IX directive 93/42/EEC)

Conformity assessment: annex VII + V (directive 93/42/EWG)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-18-03-04

Product description:

made of non-woven fabric (polyester or viscose and polyester)
with polyacrylate or hotmelt adhesive,
permeable to air and water vapour,
absorbent, non-adherent wound-pad made of viscose and polyester

Field of application / indication:

for covering minor injuries (e. g. punctures due to injections), only for external use

Product image (example):



Product variants:

Roll with skin-coloured strips (polyester), polyacrylate adhesive,
easy to tear in a hygienic dispenser.

Dispenser:

SIZE	COLOUR	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
2 x 4 cm	skin-colour	726124	04793629	4260348613539	250	132

White strips (viscose and polyester), hotmelt adhesive, loose in a box:

SIZE	COLOUR	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
2 x 4 cm	white	726127	10233225	4260348613546	500	100
2 x 6 cm	white	726128	12393068	4260348614048	300	100

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised by EO using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-18-03-08

Product description:

skin-coloured, textile fabric made of polyester with zinc oxide hotmelt adhesive, elastic, permeable to air and water vapour, absorbent rayon wound-pad made of viscose and polyester, the extended adhesive surface on one side allows the plaster to wrap around a finger or toe several times when applied and thus adhere optimally

Field of application / indication:

for covering minor injuries, only for external use

Product image (example):



Product variants:

SIZE	COLOUR	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
2 x 12 cm	skin-colour	940402	14130551	4260348616929	100	100

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-18-03-05

Product description:

white non-woven (viscose and polyester) with polyacrylate adhesive,
permeable to air and water vapour,
adapts to body contours,
transversely elastic

Field of application / indication:

for the fixation of wound dressings, cannulas or similar medical devices, only for external use

Product image (example):



Product variants:

WIDTH	LENGTH	REF	PZN (PCS)	GTIN	PCS/BOX	BOXES/CARTON
5 cm	10 m	16716001	08822670	4260348613553	1	80
10 cm	10 m	16716002	08822687	4260348613560	1	60
15 cm	2 m	16716008	14367307	4260604532321	1	128
	10 m	16716003	08822693	4260348613577	1	40
20 cm	10 m	16716004	08822701	4260348613584	1	20

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-18-02-01


Product description:

transparent, semi-permeable PU film (polyurethane) with skin-friendly acrylate adhesive,
permeable to air and water vapour,
adaptable to all parts of the body,
waterproof and watertight,
quick and hygienic to apply

Field of application / indication:

for the fixation of wound dressings, cannulas or similar medical devices, only for external use

Product image (example):



Product variants:

SIZE	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
10 cm x 10 m	719101	11860622	4260348610064	1	50
15 cm x 10 m	719102	11860639	4260348610071	1	25

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-18-02-90

Product description:

skin-coloured or white non-woven (viscose and polyester) with hotmelt adhesive, elastic, permeable to air and water vapour, absorbent, non-adherent rayon wound-pad made of viscose and polyester

Field of application / indication:

for covering minor injuries, only for external use

Product image (example):



Product variants:

SIZE	COLOUR	REF	PZN (BOX)	GTIN	ROLLS/ BOX	BOXES/ CARTON
4 cm x 5 m	skin-colour	16736010	11859530	4260348613652	1	100
	white	1673604	11859576	4260348613669	1	100
6 cm x 5 m	skin-colour	16736011	11859547	4260348613676	1	100
	white	1673605	11859582	4260348613683	1	100
8 cm x 5 m	skin-colour	16736012	11859553	4260348613690	1	96
	white	1673606	11859599	4260348613706	1	96

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-18-03-02

Product description:

skin-coloured cotton fabric with zinc oxide hotmelt adhesive,
rigid, permeable to air and water vapour,
absorbent, non-adherent rayon and polyester wound-pad

Field of application / indication:

for covering minor injuries, only for external use

Product image (example):



Product variants:

SIZE	COLOUR	REF	PZN (BOX)	GTIN	ROLLS/ BOX	BOXES/ CARTON
4 cm x 5 m	skin-colour	16738010	11859607	4260348613744	1	100
6 cm x 5 m	skin-colour	16738011	11859613	4260348613751	1	100
8 cm x 5 m	skin-colour	16738012	11859636	4260348613768	1	96

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-18-03-01

Product data sheet

Product description:

white medical tape made of artificial silk,
 with skin-friendly polyacrylate adhesive,
 tearable

Field of application / indication:

for the fixation of wound dressings, cannulas or similar medical devices, only for external use

Product image (example):



Product variants:

with side discs:

WIDTH	LENGTH	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
1,25 cm	5 m	16714003	10821471	4260348613867	12	40
	9 m	1671860	06345059	4260348613782	1	variable
		1671861	00516850	4260348613799	6	variable
		1671807	11224369	4260348613874	12	40
2,5 cm	5 m	16714004	10821488	4260348613881	12	30
	9 m	1671863	06349577	4260348613829	1	variable
		1671862	00515856	4260348613836	6	variable
		1671808	11224375	4260348613898	12	30

Product data sheet

WIDTH	LENGTH	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
5 cm	5 m	16714006	11361014	4260348613904	6	30
	9 m	16714010	13827936	4260348616783	6	30
1,25 cm	9 m	1671801	04793486	4260348613805	36	30
2,5 cm	9 m	1671802	04793575	4260348613843	24	30

without side discs:

WIDTH	LENGTH	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
1,25 cm	9 m	16714007	10180138	4260348613911	12	40
2,5 cm		1671809	10180144	4260348613928	12	30
5 cm		16714008	10180150	4260348613935	6	30

Warnings / contraindications:

- The product is not to be used in the case of an existing allergy to the material.
- The product has to be sterilised using a validated procedure before coming into contact with open wounds.
- The product has not to be used, if the packaging is damaged.
- The product has not to be reused.
- The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

- Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)
- Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)
- The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-18-01-03

Product description:

skin-coloured, textile medical tape made of cotton fabric with zinc oxide hot-melt adhesive, for normal skin, permeable to air and water vapour

Field of application / indication:

for the fixation of wound dressings, cannulas or similar medical devices, only for external use

Product image (example):

Product variants:

WIDTH	LENGTH	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
1,25 cm	5 m	16711003	10821494	4260348614017	12	40
	9 m	1671761	00516927	4260348613959	6	variable
		16711006	11224381	4260348613416	12	40
2,5 cm	5 m	16711004	10821502	4260348614024	12	30
	9 m	1671762	00517074	4260348613980	6	variable
		16711007	11224398	4260348613775	12	30
5 cm	9 m	16711008	13827942	4260348616806	6	30

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-18-01-03

Product data sheet

Product description:

white non-woven medical tape with skin-friendly hotmelt adhesive,
 for sensitive skin,
 tearable,
 permeable to air and water vapour

Field of application / indication:

for the fixation of wound dressings, cannulas or similar medical devices, only for external use

Product image (example):



Product variants:

WIDTH	LENGTH	REF	PZN (BOX)	GTIN	PCS/BOX	BOXES/CARTON
1,25 cm	5 m	1672006	11859501	4260348614055	12	40
	9 m	1672008	11224346	4260348610330	12	40
2,5 cm	5 m	1672007	11859518	4260348614062	12	30
	9 m	1672062	00517111	4260348614031	6	variable
		1672009	11224352	4260348610460	12	30
5 cm	9 m	1672010	13827959	4260348616790	6	30

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-18-01-04

Product description:

transparent medical tape made of porous polyethylene film with hotmelt adhesive,
 permeable to air and water vapour,
 can be torn lengthwise and crosswise

Field of application / indication:

for the fixation of wound dressings, cannulas or similar medical devices, only for external use

Product image (example):

Product variants:

WIDTH	LENGTH	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
1,25 cm	9 m	1671906	12393074	4260348613812	12	40
2,5 cm	9 m	1671962	00517097	4260348614079	6	variable
		1671905	11224412	4260348614093	12	30
5 cm	9 m	1671907	13827965	4260348616813	6	30

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-18-01-01

Product data sheet

Product description:

unbleached pulp from 100 % recycled material,
highly bleached pulp from TCF (= chlorine-free bleached) material,
according to DAB

Field of application / indication:

wide range of applications for maintaining hygiene standards, e.g. as protection for work surfaces, for wiping up liquids, etc., only for external use

Product image (example):



Product variants:

SIZE	COLOUR	REF	PZN (BAG)	GTIN	KG/ BAG	BAGS/ KTN
20 x 30 cm	unbleached	51123	08822730	4260348614109	5	3
	high-bleached	51423	08822753	4260348614116	5	3
40 x 60 cm	unbleached	51146	08822724	4260348614123	5	3
	high-bleached	51446	08822747	4260348614130	5	3

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

The product does not contain any substances according to REACH.

eCI@ss:


34-21-09-07

Product data sheet
Product description:

surface made of skin-friendly and hard-wearing non-woven fabric,
 highly absorbent pad made of cellulose layers,
 non-slip underlay made of PE film, which effectively prevents liquid from escaping

Field of application / indication:

suitable as incontinence underlay or as a general bed cover, or sucking up body fluids, only for external use

Product image (example):

Product variants:

SIZE	PLY	ABSORBING CAPACITY	REF	PZN (CARTON)	GTIN	PCS/ CARTON
40 x 60 cm	6	approx. 242 ml	19N03102	00527492	4260348614154	200
	12	approx. 440 ml	19N03220	00527500	4260348614161	100
60 x 60 cm	6	approx. 374 ml	19N03501	00527517	4260348614178	150
60 x 90 cm	6	approx. 572 ml	19N04027	11859719	4260348614185	50
	6	approx. 572 ml	19N04026	00527523	4260348614192	100
	12	approx. 1.122 ml	19N04037	12535930	4260348617278	50
		approx. 1.122 ml	19N04025	00527581	4260348614208	100
	20	approx. 1.837 ml	19N04024	00527598	4260348614215	50

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-39-01-90

Product data sheet

Product description:

combination of absorbent cellulose layers and a liquid-impermeable layer of PE film,
1 layer of tissue + PE film with side barrier edges

Field of application / indication:

for covering mattresses, couches, stretchers, etc. as safe protection against contamination or as a protective underlay on work tables, when taking blood samples, etc., only for external use

Product image (example):



Product variants:

SIZE	REF	PZN	GTIN	PCS/ CARTON
32 x 45 cm	1455999016	18442684	4260348614222	500
40 x 50 cm	2339028	11859808	4260348614239	1.400
80 x 90 cm	2339010	11859754	4260348614246	200
80 x 140 cm	2339013	11859760	4260348614253	100
80 x 195 cm	2339014	11859777	4260348614260	130
80 x 210 cm	2339020	11859783	4260348614277	100

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-39-01-90

Product data sheet

Product description:

combination of absorbent cellulose layers and a liquid-impermeable layer of PE film, longitudinally processed polyester threads ensure appropriate tear resistance, 2 layers of tissue + PE film with side barrier edges

Field of application / indication:

for covering mattresses, couches, stretchers, etc. as safe contamination protection or as a protective underlay on work tables, etc. if a high tear resistance and durability are required during use, only for external use

Product image (example):



Product variants:

SIZE	REF	PZN	GTIN	PCS/ CARTON
80 x 175 cm	1455999009	11859487	4260348614284	100
80 x 210 cm	23158021	11859748	4260348614291	100

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-39-01-90

Product data sheet

Product description:

surface made of skin-friendly and hard-wearing non-woven fabric,
highly absorbent pad made of cellulose layers,
non-slip underlay made of PE film that prevents liquid from escaping

Field of application / indication:

suitable as an arm underlay during dialysis, only for external use

Product image (example):



Product variants:

SIZE	PLY	REF	PZN (CARTON)	GTIN	PCS/ CARTON
22 x 55 cm	6	1850021	11859694	4260348614307	250

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.
The product has to be sterilised using a validated procedure before coming into contact with open wounds.
The product has not to be used, if the packaging is damaged.
The product has not to be reused.
The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)
Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)
The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-39-01-90

Product data sheet

Product description:

bags made of environmentally friendly paper, self-standing, space-saving storage when folded, bundled by 50 pieces

Field of application / indication:

suitable for various applications, for example as a container for storing patient-specific utensils or as a waste container

Product image (example):



Product variants:

BASE AREA	HIGHT	REF	PZN (CARTON)	GTIN	PCS/ CARTON
220 x 110 mm	285 mm	2520699	-	4260348614321	350

Warnings / contraindications:

The product has to be sterilised using a validated procedure before coming into contact with open wounds.
The product has not to be used, if the packaging is damaged.
The product has not to be reused.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

The product does not contain any substances according to REACH.

eCI@ss:

20-03-02-02

Product description:

mask according to EN 14683 made of a 3-layer,
100 % glass fibre-free non-woven fabric,
filled with a glass fibre-free filter medium,
nose clip,
elastic band or for tying,
latex-free,
virtually unhindered breathing with a bacterial filtration efficiency (BFE) > 98 %,
lightweight,
high wearing comfort thanks to anatomically correct fit

Field of application / indication:

for filtration of germs, only for external use

Product image (example):



Product variants:

SIZE	TYPE	REF	PZN (BOX)	GTIN	PCS/ BOX	BOXES/ CARTON
175 x 95 mm	TYPE II with elastic band	21221	11859725	4260348614338	50	40
	TYPE IIR to tie	21232	-	4260604534981	50	40

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has to be sterilised using a validated procedure before coming into contact with open wounds.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Notes for sterilisation:

The product may be sterilised using a validated procedure.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

Medical device of the class / rule: I / 1 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex II+III (regulation (EU) 2017/745)

The product does not contain any substances according to REACH.

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

eCI@ss:

34-15-01-02

Product description:

We are able to customise sets to meet your specific needs and requirements.

These sets can contain dressing materials, covering systems, surgical gowns, masks, bonnets or various disposable instruments.

The use of customised surgical sets offers you a number of advantages:

- Your operations are quicker and less disruptive, as the various peeling processes and jumper activities are greatly reduced.
- Your preparation time for an operation is significantly minimised thanks to the considerably simplified set-up process.
- Your stock-keeping is simplified as the number of items in stock can be reduced.
- This reduction in the number of items also relieves the workload on your logistics department, which has fewer suppliers to contact, fewer price negotiations to conduct and fewer orders and associated logistics processes to carry out.

Field of application / indication:

for various applications in the stationary or ambulatory area (see 'product variants'), for external use

Product variants:

STERILE sealed (simply packaged for stationary or ambulatory use):

EXAMPLE SETS:**Single sterile dressing set**

5 x gauze balls, 20x20cm, plum-sized, 20-thread gauze according DIN EN 14079
3 x gauze swabs, 7,5x7,5cm 8-ply, 17-thread gauze according DIN EN 14079
2 x gauze swabs, 10x10cm 8-ply, 17-thread gauze according DIN EN 14079
single packed, steam sterilised

Dressing set

1 x medical tweezer, 12,5 cm long made of plastic
4 x gauze swabs, 5x5 cm, 8-ply, 17-thread gauze according DIN EN 14079
packed in: 1 x protective pad for set packaging 60x60 cm, blue, EO-sterilised

PEG dressing set

1 x fixation plaster, 12x15cm
4 x gauze swabs, 7,5x7,5cm 8-ply
2 x gauze swabs, 7,5x7,5cm 8-ply with Y-cut
1 x instructions for use
simply packed, EO-sterilised

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

The product does not contain any substances according to REACH.

Medical device of the class / rule: Is / 4 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex VII+V (regulation (EU) 2017/745)

bzw.

Sterile Behandlungseinheit gem. Art. 12 Absatz 3 der Richtlinie 93/42/EWG

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

Product description:

We are able to customise sets to meet your specific needs and requirements.

These sets can contain dressings, draping systems, surgical gowns, masks, surgical caps or various disposable instruments.

The use of customised surgical kits offers you a number of advantages:

- Your operations are quicker and less disruptive, as the various peeling procedures and jumping activities are greatly reduced.
- Your preparation time for an operation is significantly minimised thanks to the considerably simplified set-up process.
- Your stock-keeping is simplified as the number of items in stock can be reduced.
- This reduction in the number of items also relieves the workload on your logistics organisation, which has to contact fewer suppliers, conduct fewer price negotiations and place fewer orders and associated logistics processes.

Field of application / indication:

for various applications in the surgical field (see 'product variants')

Product variants:

STERILE sealed (double-packed with counting cards for the surgical area):

EXAMPLE KITS:**Surgery kit:**

15 x laparotomy sponges, 40x40cm, 4-ply, 20-thread according DIN EN 14079, green, without loop, with X-ray
4 x laparotomy sponges, 8x80cm, 4-ply, 20-thread according DIN EN 14079, green, without loop, with X-ray
10 x gauze balls, 30x40cm, 20-thread according DIN EN 14079, with X-ray
10 x gauze balls, 25x25cm, 20-thread according DIN EN 14079, with X-ray
10 x preparation balls, 6x6cm, 24-thread according DIN EN 14079, with X-ray, in sliding box
double packed with counting card in sliding box, steam sterilised

Surgery kit urology

15 x gauze swabs, 10x13cm, 12-ply, 17-thread according to DIN EN 14079, with X-ray
20 x gauze balls, 30x40cm, 20-thread according DIN EN 14079, with X-ray, in PE bag
25 x gauze balls, 20x20cm, 20-thread according DIN EN 14079, with X-ray, in PE bag
5 x laparotomy sponges, 45x45cm, 4-fold, 20-thread according DIN EN 14079, green, with loop, pre-washed, with X-ray
10 x preparation balls, 12x12cm, 24-thread according DIN EN 14079, with X-ray, in PE bag
double-packed (wrapped in steri crepe) with counting card, EO-sterilised

Heart surgery kit

10 x laparotomy sponges, 40x40cm, 4-ply, 20-thread according DIN EN 14079, green, without loop, pre-washed, with X-ray
30 x gauze swabs, 10x10cm, 32-ply, 17-thread according DIN EN 14079, with X-ray
10 x gauze balls, 25x25cm, 20-thread according DIN EN 14079, with X-ray, in PE bag
3 x preparation balls, 6x6cm, 24-thread according DIN EN 14079, with X-ray, in sliding box
double-packed (wrapped in steri crepe) with counting card, EO-sterilised

Emergency set

20 x gauze swabs, 10x10cm, 16-ply, with X-ray
10 x laparotomy sponges, 40x40cm, 4-ply, green, without loop, with X-ray
10 x gauze balls, 30x30cm, extra-large, with X-ray
10 x gauze balls, 20x20cm, plum-sized, with X-ray
10 x preparation balls, 6x6cm, small, with X-ray, in sliding box
double packed with counting card, steam sterilised

Warnings / contraindications:

The product is not to be used in the case of an existing allergy to the material.

The product has not to be used, if the packaging is damaged.

The product has not to be reused.

The product has not to be used after the expiry date.

Storage instructions:

The product has to be protected from sunlight and stored in a dry place.

Normative and regulatory requirements:

The product does not contain any substances according to REACH.

Medical device of the class / rule: IIa / 7 (annex VIII regulation (EU) 2017/745)

Conformity assessment: article 52 (7), annex VII+V (regulation (EU) 2017/745)

or

Sterile procedure pack according to article 22 (3) regulation (EU) 2017/745

Labelling:

Labelling in accordance with DIN EN ISO 15223-1 and DIN EN ISO 20417 and with CE mark on all packaging levels.

Product data sheet

Product description:

- bag (made of 100 % cotton) with 6 different, particularly user-friendly modules with the most important First Aid items for the case of child accidents (**50 items, 83 pieces**)
- enables fast and correct action in an emergency
 - due to a structure, order and flexibility to individually identify and selectable modules sorted by topic with colour scheme
 - immediate visibility of the items contained in the modules using transparent zip bags
 - quick and easy opening of the modules using a zip lock
 - comfortable and quick opening of the bag with an non-slip, robust carrying loop
- all medical devices are hygienically packaged in zip bags and protected from contamination
- all components are labelled and described in German and English
- important information on content and areas of application as well as recommendations for action for the most common child emergencies are included in each module (module guide) as well as alphabetically sorted across all topics in the guide: First Aid for children
- with practical carrying loop for hooking onto the zip using a snap hook
- developed by pediatricians from the University Hospital of Bonn

Field of application / indication:

First Aid for child accidents, wound care for children

Product image (example):

Carrying loop is placed inside the bag; colour of the carrying loop may vary



Product variants:

SIZE OF THE BAG	REF	PZN	GTIN	BAGS/ CARTON
28 x 13 x 18 cm	452100	19105724	4260604539214	12

Content:

6 modules (50 items, 83 pieces)




WOUND CARE

Module for the treatment of minor wounds or initial treatment of larger cuts, abrasions, and lacerations

21 ITEMS (46 PIECES):

- Swabs sterile (various sizes)
- Gauze balls sterile
- Wound dressings sterile (various sizes)
- Plaster bandage
- Wound plaster for self cutting
- Adhesive tape sensitiv
- Plaster set with various plaster types
- Physiological saline solution sterile
- Scissors
- Module guide: Wound care




HYGIENE

Module to ensure hygienic First Aid measures

7 ITEMS (12 PIECES):

- Disposable gloves
- Face masks
- Multifunctional cloth
- Disinfection wipes
- Sickness bag
- Tooth rescue set Dr. Till (temporary solution for a maximum of 30 minutes)
- Module guide: Hygiene



TEMPERATURE

Module for local cold and heat applications as well as for measuring and maintaining body temperature

5 ITEMS (5 PIECES):

- Hot/cold pack
- Instant cold compress
- Rescue blanket
- Digital clinical thermometer
- Module guide: Temperature



TICKS & SPLINTERS

Module for the hygienic removal of ticks and splinters

5 ITEMS (7 PIECES):

- Ticks and splinter tweezers
- Plaster strips sensitive
- Disinfection wipes
- Magnifying glass
- Module guide: Ticks & splinters



BANDAGES

Module for the application of bandages for fixation, support, relief or compression and pressure dressings

5 ITEMS (6 PIECES):

- Conforming bandage
- Universal bandage elastic
- Triangular bandage
- Swabs sterile
- Module guide: Bandages



INFORMATION & SUPPLEMENTS

Module with helpful information material on emergencies and items for removal as well as space for individual additions

7 ITEMS (7 PIECES):

- Guidebook: First Aid for children
- Kids emergency poster
- Kids emergency card
- Kids emergency pass
- Sticker sheet
- Carrying loop (colour can vary)
- Module guide: Information & supplements

Product data sheet

Product description:

- paper bag with disposable consumables for the Dr. Till Kids Emergency Box MODULAR as a refill kit **(32 items, 65 pieces)**
- these materials should be sorted directly into the corresponding, already existing modules
- the allocation of the items to the modules can be found on the label of the refill set and on the instructions included in the refill set
- the ability to act quickly and efficiently in an emergency can only be maintained if the Dr. Till Kids Emergency Box MODULAR remains organised
- the previous versions of the Dr. Till Kids Emergency Box can also be completed with this refill set
- developed by pediatricians from the University Hospital of Bonn

Field of application / indication:

First Aid for child accidents, wound care for children

Product image (example):



Product variants:

SIZE OF THE BAG	REF	PZN	GTIN	BAGS/CTN
23 x 19 x 7,5 cm	452110	19106652	4260604539252	10

Content:


WOUND CARE

19 ITEMS (44 PIECES):

- Swabs sterile (various sizes)
- Gauze balls sterile
- Wound dressings sterile (various sizes)
- Plaster bandage
- Wound plaster for self cutting
- Adhesive tape sensitiv
- Plaster set with various plaster types
- Physiological saline solution sterile


HYGIENE

5 ITEMS (10 PIECES):

- Disposable gloves
- Face masks
- Disinfection wipes
- Sickness bag
- Tooth rescue set Dr. Till (temporary solution for a maximum of 30 minutes)


TEMPERATURE

1 ITEMS (1 PIECE):

- Instant cold compress


TICKS & SPLINTERS

2 ITEMS (4 PIECES):

- Plaster strips sensitive
- Disinfection wipes


BANDAGES

3 ITEMS (4 PIECES):

- Conforming bandage
- Universal bandage elastic
- Swabs sterile


INFORMATION & SUPPLEMENTS

1 ITEMS (1 PIECE):

- Sticker sheet

& empty ziplock bag (as a replacement)

Product description:

- mini hygiene bag made of black polyethylene with pressure closure
- tear-resistant, opaque, moisture-proof and smell-proof, can be rolled up small
- can be securely closed twice
- with adhesive strip on the back of the bag to fix the rolled-up bag
- bag size: 80 x 110 x 0,05 mm, material thickness: 50 my

available in 2 pack sizes:

- for at home in a practical sliding box with 36 bags
- for on the go in a small envelope with 6 bags

Field of application / indication:

particularly suitable for the simple, discreet, hygienic and safe disposal of tampons of all sizes, but also for the disposal of other, smaller hygiene articles (e.g. condoms, panty liners);

for at home and on the go

Product image (example):



Product variants:

SIZE	PACKAGING	GTIN	REF	PZN (PKG)	PCS/ BOX	BOXES / CTN
for at home	sliding box	4260348610019	800011	10744138	36	10
for on the go	envelope/ folding box	4260348610033	800021	10744121	6	50